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Warranty

Instead of the usual 2-year guarantee, we guarantee this device for 5 years from the date of delivery from our factory (this excludes consumables such as bulbs, tips and batteries). We guarantee the proper function of the unit provided it is used as intended by the manufacturer and according to the instructions for use. Any faults or defects which occur during the guarantee period will be repaired free of charge provided they are caused by faults in material, design or workmanship. In the case of a complaint of a defect in the product during the guarantee period, the buyer must prove that any such defect was present at the time the product was delivered. The usual legal guarantee and warranty do not apply to incorrect use, the use of non-original HEINE parts or accessories. (Especially bulbs, as these are designed with emphasis on the following criteria: Color temperature, life expectancy, safety, optical quality and performance.) Also excluded are repairs or modifications made by persons not authorised by HEINE or cases where the customer does not follow the instructions supplied with the product. Further claims, in particular claims for damage not directly related to the HEINE product, are excluded.

For **repairs**, please contact your supplier, who will return the goods to us or our authorized agent.

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Instructions for Use

HEINE DELTA 20® Dermatoscope

Read these instructions carefully before using the HEINE DELTA 20® and keep them in a safe place for future reference.

In accordance with local regulations this product should be disposed of as an electronic device separately.

Safety Instructions

Range of Applications: The HEINE dermatoscope is designed solely for the clinical examination of the surface of the skin. HEINE photo adaptor, and digital camera are specifically designed and optimized for the quick and systematic documentation of pigmented lesions and do not replace the diagnosis of the examiner.

General Safety: For documentation with a digital camera use only the HEINE photo adaptor, and our recommended adaptor from the digital camera adaptor range. Do not look directly into the light source to avoid dazzle from the intense, white, light emitting diodes (LED's). The HEINE DELTA 20® Dermatoscope is not suitable for Eye Examinations. If you use a digital camera, store your pictures first on the removable flash card and transfer the photographic images by connecting the flash card to the PC. A electrical conducting connection between camera, PC and a mains power source is not permissible.

Safety in use: Use only power sources which are designed for medical applications. Electrical and mechanical compatibility of the HEINE DELTA 20® dermatoscope can be assured only when used in conjunction with 2.5V and 3.5V HEINE battery or rechargeable handles with HEINE AV or TL-connectors.

Operation without a camera

Connect the HEINE DELTA 20® Dermatoscope to the handle. Increased mobility and ease of handling can be achieved by using the connection cable (see accessories) which connects the instrument to the handle. Prepare the skin by moistening with Dermatoscopy-oil (use a cotton swab) or disinfectant spray.

- Turn the Dermatoscope on, and place it gently so that the lesion is in the center of the contact plate (1a).
- The examiner's eye should be as close to the eye-piece (3) as possible. Adjust the focus ring until a crisp, clearly-focused image is obtained.

- For dermatoscopy, use the instrument only with the HEINE contact plate (1a, 1b) correctly-fitted.
- For general dermatological examinations, you can use the Dermatoscope without the contact plate as an illuminated loupe.
- **Tip:** Once crisp, clear focus has been achieved by adjusting the Focus ring, make a note of the position of the Focus ring (3) in relation to the linear mark (4) or edge of the instrument housing. This will help you to find your "personal" position and facilitate quicker focusing the next time you use the Dermatoscope.

Operation for digital photography

- Attach the SLR photo adaptor to the digital camera. If necessary, take a clinical photo of the lesion and general area.
- Now Adjust the Focus ring (3), aligning its lower edge with the linear mark (4). This will position the camera for the highest contrast images. The distance ring supplied with the instrument simplifies this procedure.
- To mount the photo adaptor (6) to the HEINE DELTA 20® Dermatoscope: Slide the Connector Ring (7) of the photo adaptor with one hand in the direction of the camera and hold in this position. Now guide the Dermatoscope into the photo-adaptor with the free hand until it clicks into the notch of the Dermatoscope body (5). You can now release the Connector Ring (7), it will automatically provide a tight connection between the camera and the Dermatoscope.
- Make sure that the Dermatoscope is securely connected and that it can be easily rotated inside the photo adaptor. You are now ready to start producing digital images.
- To detach the Dermatoscope from the photo adaptor slide the Connector Ring (7) on the photo adaptor in direction of the camera and pull the Dermatoscope out and away from the camera.
- **Tip:** We can only guarantee the performance of the HEINE DELTA® Dermatoscope with HEINE Photo adaptor when used with one of the following SLR digital cameras:
(www.heine.com -> DELTA 20 -> Camera for DELTA 20).
We recommend a normal visual examination of the skin before documentation with a digital camera. Please read the "Compendium of Dermatoscopy of pigmented skin lesions" supplied with the DELTA 20, which provides important information on the basic procedure.
- The resolution of a PC monitor is not high enough to display all the pixels of a modern digital camera. To get the highest resolution dermatoscopic image possible, you should increase the image size.

- The size of the field of view depends on the size of the image sensor in the camera. The field of view can vary from one digital camera to another. (www.heine.com → DELTA 20 → Camera for DELTA 20)

Brightness Control

- The voltage-regulating electronics of the HEINE DELTA 20® Dermatoscope guarantee constant brightness.
- Pressing one of the two buttons on the instrument (2) will reduce the brightness by 50% and will turn off 3 of the 6 LED's resulting in lateral illumination for improved contrast when viewing the pigmented structure.
- If the image is not bright enough, select a sensitivity of ISO 400 on the camera.

Contact Plates

- The contact plate (1a) is attached by a bayonet connection. To remove the contact plate turn it counterclockwise and pull it away from the Dermatoscope. To attach, reverse the procedure.
- For examination of hard to reach lesions use the small contact plate (1b) in place of the standard contact plate (1a). (See accessories)
- Both the standard contact plate(1a) and the small contact plate (1b) can be used for direct visualization and digital photography.

Cleaning /disinfection and sterilization

The plastics used in the instrument do not react with the dermatoscopy oil. If skin disinfectants are used, the risk of changes to the plastic, such as discolouration, cannot be completely ruled out. It is recommended that residues of this disinfectant be removed from the instrument immediately after examination in order to minimize the contact time.

Manual cleaning

- Wipe down the DELTA 20 and the photo adaptor with a damp cloth (soaked in plain water, no added agents).
It must not be rinsed in running water or immersed in water.
- Clean the lens with a lens cleaning cloth.

Cleaning and disinfecting

- Disinfectants must have been approved for use with medical equipment made of plastic in order to prevent material changes.
- Wipe down the surfaces of the DELTA 20 and photo adaptor with disinfectant.
- If you have used a spray device to disinfect the instrument, remove any excess agent from the surface in order to prevent the liquid from entering the instrument. Spray-type disinfection is only permitted with the contact plate in place.
- The internal light guide must not be treated with solutions containing alcohol.

- The contact plate can be immersed in disinfectant following removal from the instrument.

Important: The instrument must not be immersed in disinfectant.

Sterilization

Sterilize the contact plate once it has been removed from the instrument (steam sterilization: 134 °C/ 5 min.)

Technical Specifications

Magnification:	10x
Focal Range (correction):	ca. +/- 6dpt
Approved Operating Temperature:	+10 °C to 40 °C
Light Source:	High-intensity, white Light Emitting Diodes (LED), Class 1 (according to 60825-1). Due to the extreme long life of LED's they should not need replacement.

Standard and Small Contact Plates are manufactured from multi-coated Silica Glass.

Accessories

Contact Plate, without scale, Ø 23 mm	K-000.34.203
Contact Plate with graduated scale, Ø 23 mm	K-000.34.204
Small Contact Plate for hard to reach areas, Ø 8 mm	K-000.34.205
Dermatoscopy-Oil, 6 Pack (6x10 ml)	K-000.34.005
Connection Cable 1,3 mtr.	X-000.99.231
SLR Photo adaptor for Canon	K-000.34.185
SLR Photo adaptor for Nikon	K-000.34.186
SLR Photo adaptor for Olympus	K-000.34.187

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HEINE Optotechnik GmbH & Co. KG