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E mail : sav@carefitness.com

BUTTON FUNCTIONS

UP	To make upward adjustment to each function data or increase training resistance.
DOWN	To make downward adjustment to each function data or decrease training resistance.
MODE	To confirm all setting.
STAR/STOP	To start or stop workout.
RESET	To reset current setting and have the monitor switch to initial training mode for selection.
RECOVERY	To test heart rate recovery status.
BODY FAT	To test body fat % Press the BODY FAT button in standby mode and modify user data.



DISPLAY FUNCTIONS

TIME	Count up - No preset target, Time will count up from 00:00 to maximum 99:59 with each increment is 1 second. Count down - If training with preset Time, Time will count down from preset to 00:00. Each preset increment or decrement is 1 minute between 01:00 to 99:00.
SPEED	Displays current training speed. Maximum speed is 99.9 KM/H or ML/H.
RPM	Displays the Rotation Per Minute. Display range 0~15~999
DISTANCE	Accumulates total distance from 00:00 up to 99.99 KM or ML. The user may preset target distance data by pressing UP/DOWN button. Each incensement is 0.1KM or ML.
CALORIES	Accumulates calories consumption during training from 0 to maximum 9999 calories. (This data is a rough guide for comparison of different exercise sessions which can not be used in medical treatment.)
PULSE	User may set up target pulse from 0 - 30 to 230, and computer buzzer will beep when actual heart rate is over the target value during workout.
WATTS	Display current workout watts. Display range 0~999.



OPERATING PROCEDURE



Power on:

1. Please connect power adaptor to DC JACK, and console will power on with beep sound for 2 seconds and LCD display all segments.



2. Console will show "SELECT USER", user may press mode to enter into user selection mode. Use joggle wheel to select U1 to U4 and press MODE for confirmation. And then preset user information for SEX, AGE, HEIGHT and WEIGHT.

3. In standby mode, console will display “SELECT WORKOUT”, user may press MODE to enter into selection mode. And use turn joggle wheel to select MANUAL → PROGRAM → USER PROGRAM → H.R.C., → WATT







Workout in MANUAL mode:

In standby mode, select MANUAL and press MODE to enter.

Quick start: User may press START/STOP to start training in MANUAL, all exercise values will start counting up from zero.

After enter into MANUAL mode, user may set up TIME → DISTANCE → CALORIES → PULSE → RESISTANCE LEVEL by follow flashing windows, and press START/STOP to start workout. All values will start counting down to zero. (To clear setting, press RESET button)

In manual mode, biking animation will move forward every 3 km, and PC speed will be the same as user's speed.


Workout in PROGRAM mode:

In standby mode, select PROGRAM and press MODE to enter.

User may turn joggle wheel up or down to select preferred program from 1 to 12, and press MODE to confirm. Program profile will flashing, user may turn up or down to adjust profile's resistance level. TIME is fixed in 20:00, which is not adjustable. Press START/STOP to start workout and racing with PC. After start, TIME will start counting down; the runway animation will follow user's RPM input. User need to follow PC speed by checking below symbols to finish the entire 20 minutes training:

- ▲ : USER speed > PC RPM – user need to slow down
- : USER speed = PC RPM
- ▼ : USER speed < PC RPM – user need to speed up

When TIME count down to zero, console will beep for 8 seconds, and display racing result: PC WIN or USER WIN.


Workout in User program mode:


In standby mode, select USER PRO and press MODE to enter.

User may create his/her own preferred profile by turning UP and DOWN to set up resistance level of each row, and press MODE to confirm. User may hold on pressing MODE button to finish setting. Time is fixed in 20:00, which is not adjustable. User may press START/STOP button to start workout. After start, TIME will start counting down; the runway animation will follow user's RPM input. User need to follow PC speed by checking below symbols to finish the entire 20 minutes training:

- ▲ : USER speed > PC RPM – user need to slow down
- : USER speed = PC RPM
- ▼ : USER speed < PC RPM – user need to speed up

When TIME count down to zero, console will beep for 8 seconds, and display racing result: PC WIN or USER WIN.


Workout in H.R.C. mode:

In standby mode, select H.R.C. and press MODE to enter.

AGE 25 is flashing after enter H.R.C. mode, you may set your age by turning UP/DOWN button and press MODE. The monitor will calculate preset heart rate value automatically according to your age setting.

Screen will show heart rate percentage 55%, 75%, 90% and TARGET. You may select heart rate percentage by UP/DOWN/ENTER button for training.


If there is no HR input for 5 seconds, LCD will display "NEED H.R." until HR signal input.

Workout in WATT constant mode:

In standby mode, select WATT and press MODE to enter.

The preset watt value 120 is flashing on screen in WATT setting mode, select UP/DOWN/ENTER to set target value from 10 to 350. Pressing START button to start training.


BODY FAT


You may test and have the BODY FAT advice when the computer is in STOP mode. Please follow the steps as below:

1. Press BODY FAT button, and hold on handgrip to start body fat testing.
2. The symbol "-----" will display while testing period in 8 seconds. After 8 seconds, you will see the BODY FAT advice in percentage and BMI and the fat advice in different symbol.
3. It would be possible to display other indication if you have the following situation:
 - "E-1" - When you see this indication, it means you did not put your thumb properly on the conductor. Please try again.
 - "E-4" - When you see this indication, it means the BODY FAT advice is exceeded the available area which is fixed in the program.

After BODY FAT testing, press BODY FAT button to have the computer back to working continues from the previous data before testing. The BODY FAT advice figure and BMI will disappear.


SYMBOL		☐	+	▲	◆
SEX	FAT%	LOW	LOW/MED	MEDIUM	MED/HIGH
MALE		<13%	13%-25.8%	26%-30%	>30%
FEMALE		<23%	23%-35.8%	36%-40%	>40%

RECOVERY:

After exercising for a period of time, keep holding on handgrips and press "RECOVERY" button. All function display will stop except "TIME" starts counting down from 00:60 to 00:00. Screen will display your heart rate recovery status with the F1,F2....to F6. F1 is the best, F6 F6 is the worst. User may keep exercising to improve the heart rate recovery status. (Press the RECOVERY button again to return the main display.)

NOTE:

4. This computer require 24V, 2A adaptor.
5. When user stop pedaling for 4 minutes, computer will enter into power save mode, all setting and exercise data will stored until user start exercise again.
6. When computer act abnormal, please plug out the adaptor and plug in again.

Cardio – Training:

Before programming your training, it is essential to take account of your age, particularly for people of **more than 35 years of age**, as well as your physical condition. If you have a sedentary lifestyle without regular physical activity, **it is vital that you consult your doctor** to determine the level of intensity of your training. Once determined, **do not attempt to achieve your maximum during the first training sessions**. Be patient; your performance will quickly improve.

This enables you to develop your cardiovascular system (fitness of the heart / blood vessels). During Cardio-Training, the muscles need a supply of oxygen and nutritive substances. They also need to eliminate waste. Our heart increases its rhythm and sends a larger amount of oxygen to the body by means of the cardiovascular system.

The more often you exercise regularly and repetitively, the more your heart will develop, just like any other muscle. In daily life, this translates into the ability to bear greater physical effort and improve your intellectual faculties.

Definition of your work Zone:

The Maximum Heart Rate (MHR) = 220 - age (180 – age for sedentary people).


The Warming Up Zone and Cooling Down Zone are equal to 55% of your MHR.

The Fat Loss Zone is between 55 / 65% of your MHR.

The Aerobic Zone is between 65 / 85% of your MHR.

The Anaerobic Zone is between 85% of your MHR and your MHR.

Make a diagram of the target zones

The progressive training stages:

1) A warming up phase:

Start each training session by warming up progressively (**10-20 min, more if you are over 50**) which will put your muscles on standby and gently increase your heart rate. Gently increase your pulse to **55% of your MHR** (Maximum Heart Rate).

Keep your Heart Rate at this level for the duration of the warming up time. If you begin your training too quickly, you are in danger of increasing the risk of tendon or muscle injury.

2) A work phase:

This is the main part of the training session. Depending on your specific goals and your initial physical condition, weight loss, improvement of cardiovascular resistance (aerobic work), development of stamina (anaerobic work), it is possible to choose the appropriate target zone.

The Fat Loss Zone, 55 – 65%

The heart rate must be low and the duration of the exercise quite long. This zone is advised for people **resuming physical activity** after a long period of idleness, who are **trying to eliminate excess weight** or are **susceptible to cardiac problems**. Do not forget that for this type of training the duration of the exercise must be at least 30 minutes and the Heart Rate (HR) between 55 and 65% of the MHR. - Gradually increase **from 30 to 60 minutes per session**.

- Exercise **3 or 4 times a week**.

The Aerobic Zone, 65-85%

This zone is recommended for **those who are in good physical condition who train regularly**. By training in this zone, you will improve your physical condition but avoid overworking.

- Duration: **20 to 30 minutes per session**.
- Frequency: **At least 3 or 4 times a week**.

The Anaerobic Zone, 85% - MHR

This zone is exclusively for experienced athletes and is never recommended without medical approval and supervision. This zone is designed only for very fit people. It is used for interval training (or for short sprints) to help you to improve or measure your level of stamina.

The time spent in this zone must be brief to avoid the risk of accident.

3) A cooling down phase:

This entails a reduction of exertion (**10 to 20 min**). It enables you to return your cardiovascular system to normal (**55% of your MHR**) and prevent or reduce muscular pain such as cramp and stiffness. You must follow the training with a session of stretching exercises while your joints are still warm. This relaxes the muscles and helps to prevent stiffness. Stretch slowly and gently. You must never stretch until you reach the pain threshold. Hold each stretch for 30 to 60 seconds while breathing out.

Check your progress

As it improves, your cardiovascular system will have an **HR which will be lower** when at rest. You will take longer to reach your target zone and **it will be easier for you to maintain it**. Sessions will seem easier and easier and you will have more stamina in your daily life.

If, on the other hand, your **HR is higher than usual when at rest**, you will have **to rest or reduce the intensity of your training**.

Do not be discouraged by the initial sessions. Motivate yourself by scheduling specific training times. Work hard!

To order spare parts

To guarantee an efficient service, make sure you have the following items before you contact our After Sales Service:

The name or **reference of the product**.

The **manufacturing number** shown on the main chassis and on the packing box.

The **number of the missing or faulty part** shown on the exploded view of the Product described in this manual.

Contact us from Monday to Friday between 10h00 and 12h00 and 14h00 and 17h00 (except for public holidays).

C.A.R.E. After Sales Service
18 / 22 rue Bernard
Z.I. les Vignes 93012 BOBIGNY Cedex
Tel.: 0820 821 167 Fax: 00 33 1 48 43 51 03
Email: sav@carefitness.com

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