- | | |
|---|---|
| <p>Ⓓ Diagnose-Waage
Gebrauchsanweisung 2-16</p> <p>ⒼⒷ Diagnostic scale
Instruction for Use..... 17-31</p> <p>Ⓕ Pèse-personne impédancemètre
Mode d'emploi 32-46</p> <p>Ⓔ Báscula para diagnóstico
Instrucciones para el uso 47-61</p> | <p>Ⓘ Bilancia diagnostica
Istruzioni per l'uso.....62-76</p> <p>Ⓓ Diagnoz terazisi
Kullanma Talimatı77-91</p> <p>ⓇⓊⓈ Диагностические весы
Инструкция по применению..92-106</p> <p>Ⓟ Waga diagnostyczna
Instrukcja obsługi 107-124</p> |
|---|---|

Inhalt

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Contents

1	Getting to know your instrument .	18
2	Signs and symbols.....	18
3	Safety notes	19
4	Unit description	21
5	Initial use.....	21
6	Installing the control unit.....	22
7	Information	23
8	Settings	24
9	Operation	25
10	Evaluation of results	26
11	Replacing batteries	29
12	Cleaning and care of the unit....	29
13	Disposal	29
14	What if there are problems?	30
15	Technical specifications.....	31

Included in delivery

- Diagnostic scale
- Removable control unit with display
- Wall mount
- Wall mount kit
- 3 x 1.5 V type AAA (scale)
- 3 x 1.5 V type AAA (control unit)
- These instructions for use

Dear Customer

Thank you for choosing one of our products. Our name stands for high-quality, thoroughly tested products for applications in the areas of heat, weight, blood pressure, body temperature, pulse, gentle therapy, massage and air.

Please read these instructions for use carefully and keep them for later use, be

sure to make them accessible to other users and observe the information they contain. This diagnostic scale is designed to make a significant contribution to your health.

With kind regards

Your Beurer team

1 Getting to know your instrument

Function of the unit

This digital diagnostic scale is used to weigh you and to provide diagnostic data about your personal fitness.

The scale is equipped with an innovative cordless display that displays the measurements within a range of up to 2 m. The weight is displayed in 100-gram increments. The scale is intended for individual, household use.

The scale has the following diagnostic functions that can be used by up to 10 persons:

- Body weight measurement
- Measurement of body fat percentage
- Body water percentage

- Muscle percentage
- Bone mass
- Basic and active metabolic rate

The memory function allows you to access the two most recent measurements.

This scale also has the following functions:

- Switch between kilograms, pounds ("LB") and stone ("ST").
- Automatic shutoff function
- Battery change indicator for weak batteries
- When switched off, the scale displays the time and temperature.

2 Signs and symbols

The following symbols appear in these instructions.

Warning

Warning instruction indicating a risk of injury or damage to health.

Important

Safety note indicating possible damage to the unit/accessory.

Note

Note on important information.

3 Safety notes

Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain.

Safety notes

- The scales must not be used by persons with medical implants (e.g. heart pacemakers) Otherwise their function could be impaired.
- Do not use during pregnancy. Measurements may be inaccurate due to amniotic fluid.
- Do not stand on the outer edge to one side of the scale otherwise it may tip!
- Keep children away from packaging materials (risk of suffocation).

Notes on handling batteries

- Swallowing batteries can be extremely dangerous. Keep the batteries and scale out of the reach of small children. Should a child swallow a battery, seek medical assistance immediately.
- Replace weak batteries before they discharge completely.
- Always replace all batteries at the same time and use batteries of the same type.
- Batteries must not be recharged, taken apart, thrown into an open fire or short circuited.
- Leaking batteries may damage the device. If you do not intend to use the device for longer periods, remove the batteries from the battery compartment.
- If a battery has leaked, put on protective gloves and clean the battery compartment with a dry cloth.
- Batteries can contain toxins that are harmful to health and the environment. Always dispose of batteries in accordance with applicable legal regulations. Do not dispose of batteries with the normal household waste.

General notes

- The unit is for personal use only and is not intended for medical or commercial applications.
- Please note that technically related measuring tolerances are possible, as the scale is not calibrated for professional medical applications.
- The capacity of the scale is 150 kg (330 lb, 24 St). For weight measurement and bone mass measurement, the results are shown in 100 g increments (0.2 lb, 1/4 St).
- The measuring results of the body fat, body water and muscle percentage are shown in 0.1% increments.
- On delivery, the scale is set to "cm" and "kg". At the back of the scales is a switch to change between "inches", "pounds" and "stones".
- Place the scale on a firm level floor; a firm floor covering is required for correct measurement.
- Protect the unit against hard knocks, moisture, dust, chemicals, major temperature fluctuations and heat sources which are too close (stoves, heating radiators).
- Repairs may only be carried out by Beurer customer service or authorised dealers. Before submitting any complaint, first check the batteries and replace them if necessary.
- Should you have any questions concerning the use of our devices, please contact your dealer or customer service.
- All scales comply with EC Directive 2004/108/EC. Should you have any questions concerning the use of our devices, please contact your dealer or customer service.

Storage and maintenance

The precision of the measured values and service life of the unit depend on careful use:

Important

- The unit should be cleaned occasionally. Do not use abrasive detergents and never immerse the unit in water.
- Ensure that no liquid gets on the scale. Never dip the scale into water. Never wash it under running water.
- Do not place any objects on the scale when it is not being used.
- Protect the scale from impacts, moisture, dust, chemicals, large changes in temperature and sources of heat that are too close (oven, radiator). Do not press the buttons using excessive force or with pointed objects.
- Do not expose the scale to high temperatures or strong electromagnetic fields (e.g. mobile phones).

Batteries

Warning

- Swallowing batteries can be extremely dangerous. Keep the batteries out of the reach of small children. If a battery is swallowed, get immediate medical help.
- Do not throw batteries into fire. Explosion hazard!

Important

- Batteries must not be recharged or reactivated by other means, taken apart or short-circuited.
- Remove batteries if the unit will not be used for an extended period of time. This prevents damage that may occur from leaks.

Repairs


Important

- Never open or attempt to repair the unit yourself, as otherwise proper function is no longer guaranteed. Failure to observe this regulation shall void the warranty.
- If you need to have the unit repaired, contact customer service or an authorised dealer.

4 Unit description

Overview

- 1 Electrodes
- 2 Control unit with display
- 3 SET button
- 4 ▲ Up button
- 5 ▼ Down button
- 6 ON button


5 Initial use


Inserting batteries

Remove the batteries from the protective packaging and insert them into the scale and control unit according to the polarity indicated. See the diagram inside the battery compartment.

Control unit

Scale

If the scale fails to operate, remove the batteries completely and reinsert them.

For instructions on changing the batteries, see p. 29.

Changing the weight unit

The scale shows the values in kg. You can change the unit between kilograms, pounds ("LB") and stones ("ST") using the SET button. See page 24.

Setting up the scale

Set up the scale on a level and solid surface. A solid surface is a prerequisite for a correct measurement.

6 Installing the control unit

Setting up the control unit


You can either hold the control unit in your hand, set it on a level surface or keep it on the scale. The black plastic part with the Beurer logo is a magnetic wall holder for wall mounting.

Fold out the pedestal as shown in order to rest the control unit on a level surface.

Wall mounting

Using the mounting kit, you can mount the control unit on a wall. Use the mounting kit provided and install the wall holder as shown.

- Drill 2 holes in the wall next to each other horizontally, spaced 82 mm apart. Using a 6 mm drill, drill holes in the wall with a depth of approx. 40 mm. Insert the wall plugs.
- Turn the fastening screws into the wall plugs, until only about 4 mm of the screw heads protrude.
- Fasten the wall holder, then set the control unit on the magnetic holder.

Alternative wall mounting method

You can also fasten the wall holder to a wall using the double-sided adhesive tape. For this purpose, the wall surface must be clean and dry.

- Remove the protective film from the adhesive tape on the rear side of the wall holder.
- Press the wall holder firmly against the wall.

7 Information

The measuring principle

This scale operates according to the B.I.A. principle (bioelectric impedance analysis). This enables the measurement of physical relationships within seconds by means of an undetectable, completely harmless electric current. The body fat percentage and other physical relationships in the body can be determined by measuring the electrical resistance (impedance) and calculating constants and individual parameters such as age, height, gender and degree of activity.

Muscle tissue and water have good electrical conductivity, and therefore lower resistance. Bones and fatty tissue, on the other hand, have low conductivity as fat cells and bones conduct hardly any current as a result of their very high resistance.

Please be aware that values obtained from the diagnostic scale represent only an approximation of actual analytical medical data. Only a specialist physician can accurately determine body fat, body water, muscle percentage and bone structure using medical procedures (e.g. computed tomography).

General tips

- If possible, always weigh yourself at the same time of day (preferably in the morning), after using the toilet, on an empty stomach and unclothed, in order to achieve comparable results.
- Important for the measurement: only measure body fat when barefoot and with the soles of the feet slightly moist. Completely dry soles can result in unsatisfactory measurements due to inadequate conductivity.
- Stand still and upright during the measurement.

- Wait a few hours after unaccustomed physical activity.
- Wait approx. 15 minutes after getting up in the morning to allow the water in the body to distribute.
- Remember that only long-term trends are important. Short-term changes in weight within a few days are usually caused by a loss of fluid. Body water plays an important role towards our general wellbeing.

Limitations

When measuring body fat and other values, deviating and implausible results may occur in:

- children under approx. 10 years,
- competitive athletes and body builders,
- pregnant women,
- persons with fever, undergoing dialysis treatment or with symptoms of edema or osteoporosis,
- persons taking cardiovascular medication (affecting the heart and vascular system),
- persons taking vascular dilating or vascular constricting medication,
- persons with considerable anatomic deviations of the legs in relation to total body size (length of the legs considerably shortened or lengthened).

8 Settings

Enter your personal information before using the scale.

Setting the date, time and measuring unit

- With the unit in standby mode (temperature and date displayed), press the [SET] button. The date flashes in the display.
- Set the date, time and measuring unit (KG, LB, ST) one after the other by pressing the ▲ or ▼ buttons.
- Confirm each entry by pressing the [SET] button.

Setting user data

You must enter personal user parameters in order to measure your body fat percentage and other physical data.

The scale has 10 user memory locations which allows, for example, you and your family members to save personal settings and retrieve them later.

Degrees of activity

Selection of the degree of activity must refer to the medium and long term.

Degree of activity	Physical activity
1	None.
2	Low: A small amount of light physical effort (e.g. short walks, light garden work, gymnastic exercises).
3	Medium: Physical effort for 30 minutes at least 2 to 4 times a week.
4	High: Physical effort for 30 minutes at least 4 to 6 times a week.
5	Very high: Intensive physical effort, intensive training or hard physical work for at least one hour daily.

Press the [ON] button.

- Wait until the display shows "0.0".
- Press the [set] button to start changing the settings. The first memory position flashes on the display.
- Now you can enter the following settings:

User data	Parameters
Memory position	1 to 10
Body height	100 to 220 cm (3' 03" to 7' 03")
Age	10 to 100 years
Sex	male  , female 
Degree of activity	1 to 5

- Modifying values: Press ▲ or ▼ or press and hold for more rapid progress.
- Confirming entries: Press the [SET] button.




The scale is now ready for measurement.

The scale switches off automatically if no other actions are performed.

All parameters must be entered before the body fat and other values can be entered.

9 Operation

Weighing, carrying out diagnostics

1	Switching on the scale <ul style="list-style-type: none"> Press the [ON] button. <p>The display reads 0.0 KG. The scale is now ready for measurement.</p>	
2a	Weighing, weight only (without diagnostics): <ul style="list-style-type: none"> Now step onto the scale. An audible signal sounds. The display 0.0 KG starts flashing. <p>Stand still with your weight distributed evenly on both legs. After you have been weighed the results are displayed.</p> <p>When you step off the scale, the scale will switch off after a few seconds.</p>	
2b	Weighing with diagnostics: <ul style="list-style-type: none"> Select your user by pressing the ▲ button. To do so, press the button repeatedly if necessary. <p>After a short time, the display shows the ideal weight for your user data (age and height).</p> <ul style="list-style-type: none"> Now step barefoot onto the scale. An audible signal sounds. The display 0.0 KG starts flashing. <p>Stand still with your weight distributed evenly on both legs.</p> <p>After you have been weighed the results are displayed.</p>	  
<p>Important: Note that there may be no contact between your two feet, legs, calves or thighs. Otherwise the measurement will not be performed correctly.</p>		
<p>The following data are displayed automatically in sequence:</p> <ol style="list-style-type: none"> Weight and body fat percentage ("BF") with automatic evaluation, Figure 1. Weight and body water percentage "≈", Fig. 2. Weight and muscle percentage "☐", Figure 3. Bone mass "☐" and active metabolic rate "kcal", Figure 4. <p>Afterwards, weight and body fat percentage are displayed again until 0.0 appears. The scale is ready for additional measurements.</p>	   	

Retrieving measurements from memory

When the unit is switched off, the measurement is saved automatically. Two memory positions are available for each user.

<p>3</p>	<p>After a body fat measurement</p> <p>After the measurement is complete, the display reads 0.0.</p> <ul style="list-style-type: none"> Press the ▼ button. <p>The memory symbol  appears with the most recently stored measurement. The results for weight, body fat, body water, and muscle percentage are displayed.</p> <ul style="list-style-type: none"> Press the ▼ button again. <p>The second most recently saved measurement appears.</p> <ul style="list-style-type: none"> Press the ▼ button again. The memory symbol goes out and 0.0 appears in the display. 	 
<p>4</p>	<p>Switching off the scale</p> <p>The scale always switches off automatically after a certain time. An audible signal sounds.</p> <ul style="list-style-type: none"> To switch it off manually, press the [ON] button. <p>The time and temperature appear in the display.</p>	

10 Evaluation of results

Body fat percentage

The body fat values are shown in the display and evaluated with "top", "good", "average" or "poor". The following body fat percentages are guide values (contact your physician for further information).

Man

				
Age	low	normal	high	very high
10-14	<11	11-16	16.1-21	>21.1
15-19	<12	12-17	17.1-22	>22.1
20-29	<13	13-18	18.1-23	>23.1
30-39	<14	14-19	19.1-24	>24.1
40-49	<15	15-20	20.1-25	>25.1
50-59	<16	16-21	21.1-26	>26.1
60-69	<17	17-22	22.1-27	>27.1
70-100	<18	18-23	23.1-28	>28.1

Woman

				
Age	low	normal	high	very high
10-14	<16	16-21	21.1-26	>26.1
15-19	<17	17-22	22.1-27	>27.1
20-29	<18	18-23	23.1-28	>28.1
30-39	<19	19-24	24.1-29	>29.1
40-49	<20	20-25	25.1-30	>30.1
50-59	<21	21-26	26.1-31	>31.1
60-69	<22	22-27	27.1-32	>32.1
70-100	<23	23-28	28.1-33	>33.1

A lower value is often found in athletes. Depending on the type of sports, training intensity and physical constitution, values may result which are below the recommended values stated. It should, however, be noted that there could be a danger to health in the case of extremely low values.

Body water percentage

The body water percentage is normally within the following ranges:

Man

Age	poor	good	very good
10-100	<50	50-65	>65

Woman

Age	poor	good	very good
10-100	<45	45-60	>60

Body fat contains relatively little water. Therefore persons with a high body fat percentage have body water percentages below the recommended values. With endurance athletes, however, the recommended values could be exceeded due to low fat percentages and high muscle percentage.

This scale is unsuitable for measuring body water in order to draw medical conclusions concerning age-related water retention, for example. If necessary ask your physician. Basically, a high body water percentage should be the aim.

Muscle percentage

The muscle percentage is normally within the following ranges:

Man

Age	low	normal	high
10-14	<44	44-57	>57
15-19	<43	43-56	>56
20-29	<42	42-54	>54
30-39	<41	41-52	>52
40-49	<40	40-50	>50
50-59	<39	39-48	>48
60-69	<38	38-47	>47
70-100	<37	37-46	>46

Woman

Age	low	normal	high
10-14	<36	36-43	>43
15-19	<35	35-41	>41
20-29	<34	34-39	>39
30-39	<33	33-38	>38
40-49	<31	31-36	>36
50-59	<29	29-34	>34
60-69	<28	28-33	>33
70-100	<27	27-32	>32

Bone mass

Like the rest of our body, our bones are subject to the natural development, degeneration and ageing processes. Bone mass increases rapidly in childhood and reaches its maximum between 30 and 40 years of age. Bone mass reduces slightly with increasing age. You can reduce this degeneration somewhat with healthy nutrition (particularly calcium and vitamin D) and regular exercise. With appropriate muscle building, you can also strengthen your bone structure. Note that this scale will not show you the calcium content of your bones, but will measure the weight of all bone constituents (organic substances, inorganic substances and water). Little influence can be exerted on bone mass, but it will vary slightly within the influencing factors (weight, height, age, gender). no recognised guidelines or recommendations relating to bone mass measurement.

Attention:

Please do not confuse bone mass with bone density.

Bone density can be determined only by means of a medical examination (e.g. computer tomography, ultrasound). It is therefore not possible to draw conclusions concerning changes to the bones and bone hardness (e.g. osteoporosis) using this scale.

AMR

The active metabolic rate (AMR) is the amount of energy required daily by the body in its active state. The energy consumption of a human being rises with increasing physical activity and is measured on the diagnostic scale in relation to the degree of activity entered (1-5).

To maintain your existing weight, the amount of energy used must be reintroduced into the body in the form of food and drink. If less en-

ergy is introduced than is used over a longer period of time, your body will obtain the difference largely from the amount of fat stored and your weight will decrease. If, on the other hand, over a longer period of time more energy is introduced than the total active metabolic rate (AMR) calculated, your body will be unable to burn off the excess energy, and the excess will be stored in the body as fat and your weight will increase.

Results in relation to time

i Remember that only long-term trends are important. Short-term fluctuations in weight over a few days are usually the result of a loss of fluid.

The interpretation of the results will depend on changes in your overall weight and body fat, body water and muscle percentages, as well as on the period during which these changes take place. Sudden changes within days must be distinguished from medium term changes (over weeks) and long term changes (months).

A basic rule is that short term changes in weight almost exclusively represent changes in water content, whereas medium and long term changes may also involve the fat and muscle percentages.

- If your weight reduces over the short term, but your body fat percentage increases or remains the same, you have merely lost water – e.g. after a training

session, sauna session or a diet restricted only to rapid weight loss.

- If your weight increases over the medium term and the body fat percentage falls or stays the same, then you could have built up valuable muscle mass.
- If your weight and body fat percentage fall simultaneously then your diet is working – you are losing fat mass.
- Ideally you should support your diet with physical activity, fitness or power training. This way you can increase your muscle percentage over the medium term.
- Body fat, body water or muscle percentages should not be added (certain elements of muscle tissue also contain body water).

i Note regarding "ideal weight":

The ideal weight is calculated from the body height and the BMI (Body Mass Index). The ideal weight is displayed when weighing with diagnostics beginning at 18 years of age.

11 Replacing batteries

The control unit is fitted with a "change battery" indicator. If the scale is operating with batteries that are too weak, "Lo" appears on the display and the scale will switch off automatically. In this case, the batteries need to be replaced as soon as possible.

(3 x 1.5 V type AAA). The scale's batteries need to be replaced when the scale does not give a signal after switching on or does not transmit a cordless signal (3 x 1.5 V type AAA).

Note:

- When changing batteries, always use batteries of the same type, brand and capacity.
- Always replace all batteries simultaneously.
- Do not use rechargeable batteries.
- Use batteries that are free of heavy metals.

12 Cleaning and care of the unit

The unit should be cleaned occasionally.

Clean using a damp cloth, to which you can apply a little detergent if necessary. Do not use abrasive detergents and never immerse the unit in water.

Important

- Do not use aggressive solvents or cleaning agents!
- Never immerse the unit in water!
- Do not wash the unit in a dishwasher!

13 Disposal

Standard and rechargeable batteries should not be disposed of separately from the household waste. As a consumer, you are legally obliged to return used batteries for proper disposal. You can hand in your used batteries at public collection points in your district or sales outlets where batteries of this type are sold.

For environmental reasons, do not dispose of the scale with the batteries in the household waste at the end of its useful life. Dispose of the unit at a suitable local collection or recycling point.

Observe the local regulations for material disposal.

Dispose of the unit in accordance with EC Directive **2002/96/EC – WEEE** (Waste Electrical and Electronic Equipment).

Note:

The codes below are printed on batteries containing harmful substances:

Pb = Battery contains lead,
Cd = Battery contains cadmium,
Hg = Battery contains mercury.

If you have any questions, please contact the local authorities responsible for waste disposal.

14 What if there are problems?

If the scale detects an error when measuring, the following is displayed.

Display	Cause	Remedy
0_Ld	The maximum load capacity of 150 kg has exceeded.	Do not load the unit with over 150 kg.
Lo	The batteries are almost empty.	Replace the batteries (see page 21 and 29).
Err	The electrical resistance between the electrodes and the sole of your foot is too high (e.g. if you have tough calluses).	Repeat the measurement barefoot. If necessary, moisten the soles of your feet slightly. If necessary, remove the calluses from the soles of your feet.
	The fat percentage is outside the measurable range (less than 5% or greater than 50%).	Repeat the measurement barefoot or, if necessary, moisten the soles of your feet slightly.
	The water percentage is outside the measurable range (less than 36% or greater than 70%).	Repeat the measurement barefoot or, if necessary, moisten the soles of your feet slightly.
	The muscle and bone percentages are outside the measurable range (age and gender related).	Repeat the measurement barefoot or, if necessary, moisten the soles of your feet slightly.
No display	The batteries in the control unit are completely empty.	Replace the batteries (see page 21).
	The batteries in the control unit are not inserted correctly.	Check to ensure that the polarity is correct (see page 21)

No measurement possible

Possible errors	Remedy
The scale was not activated before you stood on it. If you step on the scale before "0.0 KG" is displayed, the scale will not work correctly.	Activate the scale correctly (wait for "0.0 KG") and repeat the measurement. Switch on the control unit first before the measurement. If necessary, wait 30 seconds, then repeat the measurement.
Cordless connection between scale and control unit is interrupted.	Briefly remove, then reinsert the batteries of the scale and the control unit; the scale sends an audible signal when you insert the batteries. Keep away from following electromagnetic sources: Power lines, television or computer sets, bike computers, mobile phones, weather stations
The control unit is too far away from the scale.	Reduce the distance between the scale and control unit; maximum distance is 2 metres.
The batteries in the scale are empty.	Replace the batteries (see page 21).

Incorrect weight measurement

Possible errors	Remedy
Scale is on carpeting.	Place the scale on a level and solid surface.
Incorrect zero point of scale.	Activate the scale until a signal sounds, wait approximately 4 seconds, repeat measurement.

15 Technical specifications

Dimensions of the scale	320 x 320 x 21 mm
Scale batteries	3 x 1.5 V type AAA
Control unit batteries	3 x 1.5 V type AAA
Cordless transmission	433 MHz
Measuring range	5 kg to 150 kg
Weight, bone mass unit	100 g
Body fat, water, muscle mass unit	0.1 %
BMR, AMR unit	1 Kcal

Subject to technical modifications.

The scale conforms to EC Directive 2004/108/EC + supplements.