- Ⓧ **Gebrauchsanleitung**  
Blutdruckmessgerät ..... (2–10)
- ⓊⓀ **Instruction for Use**  
Blood pressure monitor ..... (11–18)
- Ⓧ **Mode d'emploi**  
Tensiomètre ..... (19–26)
- Ⓧ **Instrucciones para el uso**  
Tensiómetro ..... (27–34)
- Ⓧ **Instruzioni per l'uso**  
Misuratore di pressione ..... (35–43)
- Ⓧ **Kullanma Talimatı**  
Bilgisayarlı tansiyon ölçer ..... (43–50)
- Ⓧ **Инструкция по применению**  
Прибор для измерения  
артериального давления ..... (51–59)
- Ⓧ **Instrukcja obsługi**  
Ciśnieniomierz ..... (59–68)

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## Dear Customer,

Thank you for choosing one of our products. Our name stands for high-quality, thoroughly tested products for the applications in the areas of heat, weight, blood pressure, body temperature, pulse, gentle therapy, massage and air. Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain.

Best regards,  
Your Beurer Team

## 1. Getting to know your instrument

The upper arm blood pressure monitor is used for non-invasive measurement and monitoring of adults' arterial blood pressure. You can use it to measure your blood pressure quickly and easily, storing the results and displaying the progression of readings together with the average.

A warning is issued for anyone suffering from cardiac arrhythmia. The values determined are classified and graphically evaluated according to WHO guidelines. Keep these instructions carefully for further use and also let other users have access to them.

## 2. Important information

### Advice on use

- Always measure your blood pressure at the same time of day, in order to ensure that values are comparable.
- Relax for approx. 5 minutes before each measurement.

- You should wait five minutes before measuring a second time.
- The readings that you take may only be used for information purposes – they are not a substitute for a medical examination! Discuss your results with the doctor. Never use them to make medical decisions independently (e.g. regarding medication and dosage thereof).
- There may be incorrect measurements where there is cardiovascular disease, as well as with very low blood pressure, circulatory disorders, dysrhythmia and other preexisting diseases.
- Only use the instrument on persons whose upper arm has the right circumference for the instrument.
- You can operate the blood pressure monitor with both batteries and mains. Note that it is only possible to save data when your blood pressure monitor is receiving power. As soon as the batteries wear out or you disconnect the mains adapter from the mains, the blood-pressure instrument loses the date and time.
- The automatic switch-off function switches off the blood pressure monitor in order to preserve the batteries, if no button is pressed within one minute.

### Storage and Care

- The blood pressure monitor is made up of precision electronic components. Accuracy of readings and the instrument's service life depend on careful handling.
  - You should protect the device from impact, moisture, dirt, major temperature fluctuations and direct exposure to the sun's rays.

- Never drop the device.
- Do not use near strong electromagnetic fields, i.e. keep it away from any radio systems and mobile phones.
- Only ever use the cuffs provided with the monitor or original replacement cuffs. Otherwise erroneous results will be recorded.
- Do not press any buttons until the cuff is in position.
- If the instrument is not used for any length of time, we recommend removing the batteries.

### Advice on batteries

- Batteries can be fatal if swallowed. You should therefore store the batteries and products where they are inaccessible to small children. If a battery has been swallowed, call a doctor immediately.
- Batteries should not be charged or reactivated with any other means, nor should they be taken apart, thrown in the fire or short-circuited.
- Remove the batteries from the instrument if they are worn out or if you are not going to use the instrument for any length of time. This prevents any damage as a result of leakage. Always replace all the batteries at the same time.
- Never use different types of battery, battery brands or batteries with different capacities. You should preferably use alkaline batteries.


### Repair and disposal

- Batteries do not belong in domestic refuse. Used batteries should be disposed of at the collection points provided.
- Never open the instrument. If these instructions are not heeded, the warranty will be null and void.
- Never attempt to repair the instrument or adjust it yourself. We can no longer guarantee perfect functioning if you do.

- Repairs may only be performed by Beurer Customer Service or authorized dealers. However, always check the batteries and replace them if necessary prior to making any complaint.
- The appliance should be disposed of according to Regulation 2002/96/EC-WEEE (Waste Electrical and Electronic Equipment). In case of queries, please contact the municipal authorities responsible for waste disposal in your area.


### 3. Unit description

1. Cuff tube
2. Cuff
3. Cuff connector
4. Cuff connector port (left side)
5. Memory button M
6. ON/OFF button 
7. Mains adapter port (rear)
8. Display
9. Date/time button 
10. Adjustment button 
11. User selection button 

### Icons in the display:

1. Date/Time Indicator
2. WHO classification
3. Memory Zones
4. Memory Sequence Number
5. Weak Battery Mark
6. Cardiac arrhythmia icon
7. Pulse Mark
8. Pulse Rate
9. Diastolic Pressure
10. Systolic Pressure
11. Problem shown on the display

Note: Batteries containing pollutant substances are marked as follows: Pb = Battery contains lead, Cd = Battery contains cadmium, Hg = Battery contains mercury.

### Setting date and time

It is vital to set date and time. Only in this way it is possible to save and subsequently retrieve your measured values with the right date and time.

Press the “⊖” key and “⊕” to adjust the month. Then press the “⊖” and “⊕” in turn to adjust the date, hour, minute and confirm the entry.

Time is shown in 12 hour format.

## 4. Prepare measurement

### Inserting battery

- Remove the battery cover from the back of the monitor.
- Insert four AA 1.5 V alkaline batteries, making absolutely sure that you insert them with the correct polarity as marked. Never use rechargeable batteries.
- Replace the battery cover carefully.

If the battery change „“ is continuously illuminated, measurement is no longer possible and you must replace all the batteries. Once batteries have been removed from the device, the time must be reset.

Used batteries do not belong in the household waste. You are legally obliged to dispose of the batteries. Dispose of them via your specialist electrical supplier or local collecting point for recyclable waste.


### Mains operation

You can also operate this device with a mains adapter. No batteries should be in the battery compartment for this. The mains adapter is available from retailers or from the service address under order no. 071.03. The blood pressure monitor should only be used with the mains adapters specified here. The mains adapter should only be connected to the mains voltage indicated on the rating plate. Once you have disconnected the device from the mains, the date and time on the blood pressure monitor are lost. Any measurement results that have been stored are however retained.

## 5. Measuring blood pressure

### Positioning cuff




Fit the cuff round your bare left upper arm. Blood circulation in the arm should not be restricted by tight clothing or other objects.

The cuff should be placed on the upper arm so that the lower edge is 2 to 3 cm above the bend of the elbow and above the artery. The tube should be in line with the centre of the palm.




Now place the free end of the cuff snugly, but not too tightly, around the arm, and fix it with the Velcro fastener. The cuff should be fitted tight enough to allow just two fingers to fit beneath the cuff.

Insert the cuff tubing into the socket for the cuff attachment.

**Important:** The instrument should only be operated with the original cuff. The cuff is suitable for an arm circumference of 24 to 36 cm. A larger cuff for upper arm circumferences from 34 to 46 cm is available from retailers or the service address under order no. 162.383.

### Correct posture


- Rest for approx. 5 minutes before each measurement. Otherwise there may be divergences.
- You can perform the measurement either sitting or lying down. Always make sure that the cuff is on a level with your heart.
- In order not to distort the result, it is important to keep still during the measurement and not talk.

### Select memory

You have two memories (50 memory spaces each) in order to save the test results of 2 different people separately, or else save measurements in the morning and evening separately. Select the required memory space by pressing the user selection button .

### Measuring blood pressure

- Put on the cuff as described previously and assume the position in which you want to carry out the measurement.
- Press the “User-Switching key”  to select a memory zone, and press the “⊖” key 2x to start measurement in the chosen memory zone. After checking the display with all digits lighting up, the monitor will automatically inflate. Following the self-test, during which all display elements briefly appear, the measurement begins. During inflation the device already determines measured values used to estimate the required inflation pressure. If this pressure is insufficient, the device automatically inflates by another 40 mmHg (Real Fuzzy Logic).
- Then the pressure in the cuff is slowly released and the pulse is measured.
- When the measurement has been completed, the remaining air pressure is released very quickly. The pulse, the systolic and the diastolic blood pressure are displayed.

- The device requires approx. 5 seconds to store the determined measured value.
  - You can interrupt measurement at any time by pressing the ON/OFF button .
  - To switch off and release the pressure, press the ON/OFF button  again. If you forget to switch off the device, it switches off automatically after approx. 1 minute.
- Wait at least 5 minutes before taking another measurement!

## Evaluating results

### Cardiac arrhythmia:

This instrument can identify possible cardiac arrhythmia disorders during measurement and if necessary indicates the measurement with the flashing icon .

This may be an indicator for arrhythmia. Arrhythmia is a condition where the heart rhythm is abnormal as a result of defects in the bioelectrical system controlling the heart beat. The symptoms (omitted or premature heart beats, slow or excessively fast heart rate) may be caused, among other things, by heart disease, age, physical predisposition, excessive use

of stimulants, stress or lack of sleep. Arrhythmia can only be ascertained through examination by your doctor. Repeat the measurement if the flashing icon  is displayed after the measurement. Please note that you should rest for 5 minutes between measurements and not talk or move during the measurement. If the icon  appears often, please contact your doctor. Any self-diagnosis and treatment based on the test results may be dangerous. It is vital to follow your doctor's instructions.

### WHO classification:

The World Health Organization (WHO) and National High Blood Pressure Education Program Coordinating Committee have developed a blood pressure standard, according to which areas of low- and high-risk blood pressure are identified. This standard, however, is a general guideline as individual's blood pressure varies among different people and different age groups...etc.

It is important that you consult with your physician regularly. Your physician will tell you your normal blood pressure range as well as the point at which you will be considered at risk.

**Table for classification of blood pressure values (unit of measurement mmHg) for adults:**

Range	Systolic	Diastolic	Action
Hypotonia (degraded blood pressure)	lower than 105	lower than 60	Check with doctor
Normal range	between 105 and 120	between 60 and 80	Self check
Pre-hypertensive blood pressure range* P	between 120 and 140	between 80 and 90	Check with doctor
High blood pressure - degree 1	between 140 and 160	between 90 and 100	Consult your doctor
High blood pressure - degree 2	higher than 160	higher than 100	Consult your doctor

(\* ) Blood pressure range that might transfer to high blood pressure

Adapted from JNC 2003


The WHO classification in the display shows the range of the measured blood pressure.

If the values for systolic and diastolic pressure are in two different WHO ranges (e.g. systolic in the grade 1 hypertension range and diastolic pressure in the normal range), the WHO classification on the unit always indicates the higher range (grade 1 hypertension in the example described).

## 6. Saving, retrieving and deleting results

The computer automatically saves the blood pressure readings of the last 50 measurements. When the 50 memory locations are exceeded, the current oldest value is deleted.

You can display the memory as follows:

- Press the “User-Switching key”  to select a memory zone, and press the “M” key to switch from memory to memory with sequence numbers.
- The saved measured values are continuously stored.
- Press the On/Off switch  to turn the unit off.
- If you forget to switch off the computer, it will switch off automatically after 1 minute.

Press the “User-Switching key”  to select a memory zone, then press and hold the “M” key for approximately 5-seconds to clear the data in the pre-designated memory zone. The data in both memory zones – zone 1 and zone 2 may be erased if any of the batteries is removed.

## 7. Cleaning and storing the instrument

- Clean your blood pressure computer carefully only with a slightly moistened cloth.
- Do not use detergents or solvents.

- On no account must you immerse the computer in water, otherwise liquid can enter it and cause damage.
- When storing the device, make sure that no heavy objects are placed on top of it. Remove the batteries. The cuff tube should not have any sharp kinks.

## 8. Rectifying faults

Error messages can occur when

- the blood-pressure value is unusually high or low („EE“ appears in the display),
- you move or talk during the measurement („EE“ appears in the display),
- the cuff tube is not properly inserted („E1“ appears in the display),
- inflation takes longer than 20 seconds („E1“ appears in the display),
- the inflation pressure is higher than 300 mmHg („E2“ appears in the display),
- an error occurs when storing the measured values („E3“ appears in the display).
- Exceeding measurement range („Er“ appears in the display).

In the about cases, you must repeat the measurement. Make sure that the cuff tube is properly inserted and that you do not move or talk. Put the batteries back in if necessary, or else replace them.

## 9. Specifications

Model No.	BM 16
Method of Measurement	Oscillometric, non-invasive blood pressure measurement on the upper arm

Range of Measurement	Pressure 0–299 mmHg, systolic 40–250 mmHg, diastolic 40–250 mmHg, pulse 40–199 beats/minute
Accuracy of the pressure reading	systolic $\pm 3$ mmHg / diastolic $\pm 3$ mmHg / pulse $\pm 5\%$ of the displayed value
Measuring uncertainty	max. permissible standard deviation according to clinical testing: systolic: 8 mmHg / diastolic: 8 mmHg
Sets of Memory	2 x 50 memory places
Unit Dimension	L156 x W105 x H74,5 mm
Unit Weight	Approx. 330 g
Cuff Size	Wrist circumference approx. 24–36 cm
Operating Temperature	+10 °C to +40 °C, 40-85 % RH
Storage-Temperature	-10 °C to +60 °C, 10-90 % RH
Power Supply	4 x 1,5 V AA batteries (alkaline type LR6)
Battery Life	For approx. 300 measurements, depending on blood pressure level or inflation pressure
Accessories	With storage pouch, instruction manual, 4 pcs “AA” battery
Classification	Application part Type BF
Key to symbols:	Important! Read the operating instruction! 
	Application part Type BF 

These specifications are subject to change without notice for purpose of improvement.

## 10. Adapter

Model No.	AN 5706
Input	230V, 50Hz
Output	12V DC, 800mA, only in connection with beurer blood pressure monitor
Supplier	Ansmann AG
Protection	This device is double insulated and protected against short circuit and overload by a primary thermal fuse. The green LED lights under normal conditions. Make sure to take the batteries out of the compartment before using the adapter
	double insulated / equipment class 2
Enclosures and Protective Covers	Equipment enclosed to protect against contact with live parts, and with parts which can become live (finger, pin, hook test) The operator shall not contact the patient and the output plug of AC adaptor simultaneously.

- This unit is in line with European Standard EN60601-1-2 and is subject to particular precautions with regard to electromagnetic compatibility (EMC). Please note that portable and mobile HF communication systems may interfere with this instrument. For more details, please contact customer service at the address indicated.

- The instrument is in line with the EU Medical Products Directive 93/42/EC, the German medical products act and European standards EN1060-1 (non-invasive blood pressure monitors, Part 1: general requirements) and EN1060-3 (non-invasive blood pressure monitors, Part 3: Supplementary Requirements for Electromechanical Blood Pressure Measurement Systems).
- If you are using the instrument for commercial or economic purposes, you must check measurement precision regularly in accordance with the 'Operators' Ordinance on Medical Products'. Even in the case of private use, we recommend checking measurement precision every two years at the manufacturers.

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