

**Codefree**

BLOOD SUGAR  
MEASURING DEVICE

Step by step





## Contents

<b>1 About Your Beurer Blood Sugar Measuring Device</b> .....	<b>4</b>
1.1 Delivery scope, replacement parts and accessories .....	5
1.2 Functions of the unit .....	6
1.3 Signs and symbols .....	7
<b>2 Safety information and warnings</b> .....	<b>8</b>
<b>3 Description of the unit and accessories</b> .....	<b>12</b>
3.1 Blood sugar measuring device .....	12
3.2 Lancing device and needle lancets .....	12
3.3 Display symbols .....	13
3.4 Test strips.....	14
<b>4 Using the unit for the first time and basic settings</b> .....	<b>16</b>
4.1 Removing the battery insulation strip, replacing the batteries .....	16
4.2 Basic settings .....	17
<b>5 Measuring your blood sugar</b> .....	<b>18</b>
5.1 Acquiring blood samples .....	18
5.2 Measuring the blood sugar value .....	22
5.3 Evaluating a blood sugar measured value.....	23
5.4 Function check with control solution.....	25
<b>6 The Measured Value Memory</b> .....	<b>28</b>
6.1 Displaying individual values.....	28
6.2 Displaying average blood sugar levels .....	29
6.3 Deleting the measured value memory .....	29
6.4 Transferring measured values to a computer .....	30
<b>7 Storage and maintenance of the unit</b> .....	<b>32</b>
<b>8 Troubleshooting</b> .....	<b>33</b>
<b>9 Technical specifications</b> .....	<b>35</b>
<b>10 Guarantee and Customer Service</b> .....	<b>38</b>

## 1 About Your Beurer Blood Sugar Measuring Device

Dear valued customer,

Thank you for choosing one of our products. Our name stands for high-quality, thoroughly tested products for applications in the areas of heat, weight, blood pressure, body temperature, pulse, gentle therapy, massage and air.

Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain.

With kind regards  
Your Beurer team

### About Your Beurer Blood Sugar Measuring Device

The blood sugar measuring device is designed for fast, uncomplicated measurement of blood sugar.


You can use it to determine your blood sugar content, save the measured values and display the average of the measured values.

The large display with background lighting makes the measured values easy to read. The user-friendly design with the handy measurement strips and reduction of the operating controls to 2 buttons guarantees simple, yet accurate measurements.

The unit can be connected to a computer via a USB cable (not included in the scope of delivery). Using special software, you can evaluate your measured values on the PC and use the evaluations for monitoring your blood sugar levels.

## 1.1 Delivery scope, replacement parts and accessories

Inspect your set for external damage to the carton packaging and for completeness of the content.



Item	Designation
A	1 blood sugar measuring device
B	1 lancing device with AST cap for taking blood samples from alternative sites
C	10 test strips
D	10 sterile needle lancets
E	2 button cell batteries 3 V CR2023 (pre-installed)
F	1 practical carry case
	These instructions for use, additional information material

- If the box has sustained extensive damage or if any contents are missing, please return the system to the retailer.
- The blood sugar measuring device, test strips and control solution (available to purchase separately) are all designed to be used together. Therefore use only the test strips and control solution intended for this measuring device.



### Note

- Only use the manufacturer's original accessories.

## Purchasing supplies

Test strips, control solution and lancets are also available for purchase without a prescription.

Item	REF
50 test strips	REF 464.01
LEVEL 1 control solution	REF 463.05
100 needle lancets	REF 457.01
Beurer GL40 PC Kit (connection cable, driver and test software)	REF 463.10

## 1.2 Functions of the unit

---

This unit is designed for measuring the blood sugar content. It is suitable for household use.

Using this measuring device, you can quickly and easily:

- Measure your blood sugar.
- Display and save your measured values.
- Display your average blood sugar measurements from the last 7, 14, 30 and 90 days.
- Set the time and date.
- upload your stored measured values to a PC for evaluation (additional accessories required).

The measuring device also features the following control functions:

- A warning in the event of unsuitable temperatures.
- Battery change reminder when battery power is low.















### Warning




- **This unit is only to be used for regular monitoring and not for the diagnosis of diabetes.**
- **Discuss your insulin dose with your doctor.**

### 1.3 Signs and symbols

The symbols on the packaging and the nameplate of the measuring device and accessories have the following meaning:

	In-vitro diagnostics		Manufacturer
SN:	Serial number		Refer to instructions for use
	Temperature range +2°C to +30°C		PCT: Certification for products exported to the Russian Federation and CIS countries
	Not for reuse/For single use only		Grüner Punkt: German recycling system
	Expiry date		Contents sufficient for <n> tests
	Maximum shelf life after opening in months	REF/Art.-Nr.	Order number
LOT	Lot number		Biohazard, risk of infection
STERILE R	Sterilised by radiation (lancets)	mg/dL	Unit of measure for blood sugar value
	Caution, refer to accompanying documents		

The following symbols in these instructions mean

-  **Warning**      Warning instruction indicating a risk of injury or damage to health.
-  **Important**      Safety note indicating possible damage to the unit/accessories.
-  **Note**              Important information to be noted.

## 2 Safety information and warnings



### Risk of infection

All components of the measuring device and accessories may come into contact with human blood, and are therefore a possible source of infection.

### Warning

- **This measuring device must display the blood sugar content in mg/dL. The measuring unit mg/dL appears next to each blood sugar level reading. Be sure to contact customer service if your measuring device does not display mg/dL. You risk damaging your health if you perform a blood sugar measurement using an unfamiliar unit of measure, since values may be interpreted incorrectly and cause the wrong corrective measures to be taken.**
- This measuring device must only be used by one person. There is a risk of infection if the same unit is used by more than one person.
- The lancing device is suitable for household use. Never share the lancing device or needle lancet with other people (**risk of infection!**).
- Use a new sterile needle lancet for each blood test (**for one-time use only**).

### General notes

#### Warning

- Do not use the unit near strong electromagnetic fields, wireless systems or mobile telephones.

### Measuring blood sugar

#### Warning

- The measured values you determine are useful for information purposes only – they are not intended to replace consultation with your doctor. Discuss your measured values with your doctor regularly. Do not change any aspects of your treatment unless instructed by your doctor.
- Water deficiency or a large loss of fluid, for example by sweating, can lead to false measurement results.
- A very high or very low haematocrit value (proportion of red blood cells) can lead to incorrect measurements. With a very high haematocrit value (over 60%), the displayed blood sugar value may be too low. If you have a very low haematocrit value (below 20%), the blood sugar value may be too high. If you do not know your haematocrit value, ask your doctor.
- Do not use the test strips to take blood sugar measurements of newborn infants.
- Metabolites such as ibuprofen, sodium salicylate, tetracycline, tolbutamide, unconjugated bilirubin, cholesterol, creatinine, triglycerides, galactose, maltose, xylose, paracetamol, gentisic acid, levodopa, dopamine, methyl dopa, uric acid and ascorbic acid do not affect the result so long as they lie within the physiological value range.
- Do not use this unit to test any seriously ill patient.



### **Warning**

- Lipemia effects: Increased blood triglycerides up to 3000 mg/dL hardly affect the results. Above this level, however, the blood sugar test may be affected.
- Only use fresh capillary whole blood. Do not use serum or plasma.
- Use capillary blood obtained without squeezing the puncture site. Squeezing causes interstitial fluid to dilute the blood and leads to a false measurement result.
- Do not use the test strips at elevations over 3048 m.

### **Note**

- The Beurer GL40 mg/dL measuring system is suitable for measuring capillary whole blood.

### **Storage and maintenance**

#### **Warning**

- Store the measuring device and accessories out of the reach of young children. Small parts, such as needle lancets, batteries or test strips, can be life-threatening if swallowed. If parts are swallowed, seek medical advice immediately.
- The test strip container contains a desiccant that may cause skin or eye irritations if inhaled or swallowed. Keep the container out of reach of young children.

The measuring device is made from precision and electronic components. The accuracy of the measured values and service life of the unit depend on careful use:

- Protect the unit and accessories from impacts, moisture, dirt, large temperature fluctuations and direct sunlight. Do not store the unit, the test strips and the control solution in a refrigerator, bathroom or car!
- Do not let the unit fall.

### **Batteries/Saving measured values**

#### **Warning**

- Make sure batteries are not accessible to children. Children can put batteries in their mouth and swallow them. This can cause severe harm to their health. In this case, consult a physician immediately!
- Normal batteries must not be recharged, heated or thrown into an open flame (**danger of explosion!**).

## Important

- Do not disassemble or short-circuit the batteries.
- Always replace all batteries at the same time and use batteries of the same type. Do not use any rechargeable batteries.
- Leaking batteries may damage the unit. If you do not intend to use the unit for a prolonged period, remove the batteries from the battery compartment.

### Caution!

- Leaking or damaged batteries can cause burns if they come into contact with your skin. In this case, wear suitable safety gloves.

## Note

- The stored blood sugar measured values are retained when you replace the batteries. The date and time are also retained when you replace the batteries and when the batteries are empty.
- Only use lithium-ion batteries.

## Repairs

### Note

- Never open the unit. Opening the unit invalidates the guarantee.
- Do not perform any repairs on the unit yourself, otherwise correct functioning of the unit can no longer be guaranteed.
- If repairs are required, please contact customer service.

## Disposal

### Warning

- When disposing of the materials of the measuring device, be sure to observe the generally applicable safety precautions for handling blood. Carefully dispose of all blood samples and materials with which you have come into contact, in order to prevent injury and infection of others.
- After use, dispose of the test strips and lancets in a sharp-proof container.



## Note



In the interest of protecting the environment, the unit must not be thrown out with the household waste at the end of its service life. Please dispose of the unit in accordance with EU directive 2002/96/EC – WEEE (Waste Electrical and Electronic Equipment). For any queries, please contact the municipal authority responsible for disposal.

Standard and rechargeable batteries should not be disposed of along with household waste. As a consumer, you are legally obliged to return used batteries for proper disposal. You can hand in your used batteries at public collection points in your district or sales outlets where batteries of this type are sold.

The following symbols are found on batteries that contain harmful substances:

Pb = battery contains lead


Cd = battery contains cadmium

Hg = battery contains mercury




### 3 Description of the unit and accessories

#### 3.1 Blood sugar measuring device




Item	Designation
1	Receptacle for test strips, with lighting
2	Display
3	ON/OFF button
4	Computer port
5	Battery compartment (underside)
6	▲▼ rocker button

#### 3.2 Lancing device and needle lancets



Item	Designation
1	Cap
2	Protective disk for the lancet
3	Sterile needle lancet
4	Lancet holder
5	Release button
6	Tensioning device

### 3.3 Display symbols



- |   |  |    |                       |
|---|--|----|-----------------------|
| 1 | Date   | 7  | Battery change symbol |
| 2 | Time   | 8  | Symbol to apply blood |
| 3 | Computer connection active                                   | 9  | Temperature symbol    |
| 4 | Measured value display, HI / LO display, average blood sugar | 10 | Blood drop symbol     |
| 5 | Blood sugar unit mg/dL                                       | 11 | Test strip symbol     |
| 6 | Blood sugar unit mmol/L – not functional                     | 12 | Memory symbol         |
|   |  | 13 | Error symbol          |
|   |  | 14 | Speaker symbol        |




#### Note

The measuring device is supplied with the following basic settings:

- Blood sugar unit: mg/dL
- Beep on

### 3.4 Test strips

Front view




- 1 Slit for blood collection
- 2 Hold here
- 3 Contacts

Rear view



You will know which side is the rear from the contact panels.

Insert the test strip into the unit so that the contacts are visible in the slit. Ensure that the front side of the test strip is facing you.



### **Note**






Read the following information about handling and storing your test strips carefully. The test strips are guaranteed to yield accurate measurement results only if you observe all instructions.

### **Warning**

- Each test strip is to be used only **once** and only for **one** patient!

## Handling test strips

### **Note**

- Once you have removed the test strip, tightly close the container again immediately.
- Do not use test strips if the expiry date has passed. The use of expired test strips can lead to inaccurate measurements. The expiry date is located on the container next to the hourglass symbol .
- The test strips can be kept for three months once the container is open. Make a note of when this time will run out (date of opening + 3 months ) on the label. The shelf life decreases if the expiry date is reached before the end of the 3 months (see the date next to the hourglass symbol .
- Do not use the test strips if either of these dates has expired (  /  ).
- Any part of the test strip may be handled with dry, clean hands.
- Use the test strips immediately upon removal from the container.
- Do not bend, cut, or otherwise modify the test strip.
- Do not use test strips that have come into contact with liquids for the blood sugar measurement.

## Storing test strips

### **Note**


- Store test strips in a cool, dry place between +2°C and +30°C. Protect test strips against direct sunlight or heat. Do not store test strips in the fridge.
- Permitted relative humidity under 90%.
- Store test strips only in the original container – never use another container.

## 4 Using the unit for the first time and basic settings

### 4.1 Removing the battery insulation strip, replacing the batteries


#### Note

- Two batteries are included in the delivery scope of your blood-sugar measuring device. They are already inserted in the battery compartment.
- Before using the unit the first time, you must remove the insulation strip.



- 1 Remove the cover of the battery compartment on the underside of the unit.
- 2 When you perform a battery change, remove all batteries. The unit keeps the date and time while you change the batteries, if there is at least one battery inserted. If the battery change lasts a long time, reset the date and time (see "Basic settings", page 17).
- 3 Insert two new batteries of type **CR 2032 3V**.  
Ensure that the batteries are inserted with the correct polarity.  
See the diagram inside the battery compartment.
- 4 Carefully replace the lid on the battery compartment.

#### Note

- When the battery change symbol  appears, the batteries are almost empty. Replace both batteries as soon as possible.
- When "LP" is displayed, the batteries are so empty that no more measurements are possible.



## 4.2 Basic settings

- 1 Remove the batteries and then reinsert them. A beep sounds.  
The year display flashes.



2011

- 2 **Setting the date and time.**

 **Note**

- You must set the date and time. Only if you do that can you save your measured values correctly with the date and time and later retrieve them.
- The time is displayed in 24-hour format.

Set the year (calendar extends to 2099) by pressing the ▲ or ▼ button. Confirm via the ON/OFF button [3].

The day display flashes.

Repeat the above procedure for the day, month, hours and minutes.

"0n" and the speaker symbol are displayed.

- 3 **Turning the beep on/off**

To turn off the beep, press the ▲ or ▼ button. "OFF" is displayed.

The speaker symbol disappears from the display.

Confirm via the ON/OFF button [3].

- 4 The measuring device switches off automatically.

## 5 Measuring your blood sugar

### Warning

If you drop the lancing device with a needle lancet inserted, carefully pick it up and dispose of the lancet.

### Important

- Use the lancing device only with needle lancets from the manufacturer. Use of other needle lancets can impair the function of the lancing device.
- If the lancing device is manufactured by a third party, consult their instructions for use.

### 5.1 Acquiring blood samples

---

#### Preparing to take the blood sample

- 1 Choose a site on the body from where you wish to take the blood sample. You can use the lancing device to get blood samples from the fingertip or other parts of the body, such as the palm, forearm or upper arm. We recommend that you take the blood sample from the fingertip. To make the blood sample as pain-free as possible, do not take the blood directly from the centre of the fingertip, but instead from slightly to the side of the centre.

### Warning

- **If you suspect hypoglycaemia: be sure to draw blood from the fingertip.** This is because changes in the blood sugar level can be measured quickly in blood samples from the fingertip.
  - Measuring from a fingertip and another part of the body (AST) can lead to significantly different measured values. Always consult your doctor before starting to test at alternative sites.
- 2 Have the following components ready: Measuring device, container of test strips, lancing device, and a sterile needle lancet. If obtaining a blood sample from a site other than the fingertip, you will also need the AST cap.
  - 3 Before taking the sample, wash your hands using soap and warm water. In addition to optimal hygienic conditions, this also ensures good circulation of blood through the fingertips. Dry your hands carefully. Also ensure that your lancing site is hygienically clean if taking a blood sample from an alternative site (AST).

### Warning


If you wipe the lancing site with alcohol, ensure that the site has dried completely before beginning the measurement.

## Taking the blood sample


### Warning

- Change the prick point for each test, e.g. another finger or the other hand. Repeated pricks at the same point can lead to inflammation, numbness or scars.
- Do not use the AST cap to take a blood sample from the finger.
- Never squeeze the finger to obtain a larger droplet of blood. Squeezing causes interstitial fluid to dilute the blood; this can lead to a false measurement result.
- Note that inadequate circulation at the puncture site, for example caused by cold or illness, can lead to inaccurate measurements.


1 Twist to remove the cap from the lancing device.



2 Place a sterile needle lancet in the lancing device and push the lancet firmly into place.



3 Remove the protective disc from the lancet by twisting while holding the shaft of the lancet firmly. Retain the protective disc so that you can safely dispose of the used needle lancet after taking the blood sample.




4 You will need to use a different cap depending on the site from which you want to take the sample:

**Fingertip:** Cap (white)

**Other parts of the body:** AST cap (transparent)

Place the selected cap on the lancing device and twist it into place.




5 **Setting the puncture depth:**


Five different puncture depths can be selected on the lancing device:

- 1 to 2: Delicate or thin skin
- 3: Normal skin
- 4 to 5: Thick or rough skin

Turn the cap in the corresponding direction until the arrow points to the desired puncture depth.




- 6 Pull back the tensioning device until you hear it click into place. If it does not engage, the lancing device may have unintentionally already been tensioned when inserting the needle lancet and could already be engaged.



- 7 The lancing device can now be used to take the blood sample. Ensure that the blood remains in a droplet form and is not smeared.


#### Blood sample from the fingertip

Place the lancing device firmly against the finger, slightly to the side of the centre of the fingertip. Press the release button. Lift the lancing device from the finger. You need a round blood drop of at least 0.6 microlitres to form (corresponds to approx. 1.4 mm, original size: ●).




#### Blood sample from an alternative site (AST)




Choose a site that is soft, not next to a bone, with no visible veins and without excessive body hair. Warm the puncture site to ensure good circulation, for example by lightly massaging the area. Place and hold the lancing device against the puncture site for a few seconds and press the release button. Continue to hold the lancing device against your skin until a round droplet of blood has formed under the cap. Keep applying pressure until the size of the blood drop is at least 0.6 microlitres (corresponds to approx. 1.4 mm, original size: ●). Lift the lancing device carefully away from the skin.





- 8 If insufficient blood was collected, increase the puncture depth and repeat steps 5 to 7.
- 9 Untwist the cap and carefully remove it from the lancing device.





- 10 Place the protective disk you kept earlier flat on a hard surface. Pierce the protective disk with the tip of the needle so that the needle is covered.
- 11 Carefully remove the needle lancet from the lancing device and discard the lancet in a suitable sharp-proof container. Carefully dispose of all blood samples and materials with which you have come into contact. This prevents injury and prevents infection spreading to other people.
- 12 Twist the (white) cap back in place.



## 5.2 Measuring the blood sugar value

- 1 Hold the measuring device so that the display is facing you.
- 2 Insert a test strip, contacts first, into the unit. Ensure that the front side is facing you. You can handle any part of the test strip with clean, dry hands.
- 3 The unit switches on automatically and shows the start display. As soon as the hand  and  symbols are flashing, the unit is ready to perform measurements.
- 4 Hold the blood collection slit (at the point of the test strip) against the blood droplet. Do not press the puncture site (fingertip or other body part) against the test strip. The blood must not be smeared. The blood is sucked into the slit.



-  **Warning:** Hold the blood collection slit of the test strip against the blood droplet until the slit is completely filled and you hear a beep. Withdrawing the test strip from the blood droplet before the beep may lead to incorrect measurements.
- 5 Once the slit is filled with blood, the unit performs the blood sugar measurement. The measuring device counts down for approx. five seconds. Then the measurement result is shown in the display.  
Read the measured value.  
For explanations and corrective measures concerning measured values, see the next Chapter "Evaluating a blood sugar measured value", page 23.  
If an error message is displayed, read the Chapter "Troubleshooting", page 33.
  - 6 Remove the test strip from the unit and dispose of it carefully in accordance with the applicable guidelines in order to prevent infection of other people.

### **Note**

- If the unit does not begin the measurement, do **not** subsequently apply more blood. Instead, remove the test strip and end the test procedure. Start again using a new test strip.
- If the test strip is already in the unit and you do not put any blood on it within two minutes, the unit switches off. Remove the test strip then and reinsert it into the slot, so that the unit switches on again automatically.
- If you experience difficulties applying blood to the test strip correctly, contact customer service.
- If you are measuring in a dimly lit room, press the ON/OFF button to switch on the unit. The test strip area lighting turns on and makes it easier for you to insert the test strip. In addition, the backlighting turns on for the results display.

### 5.3 Evaluating a blood sugar measured value

---

Your blood sugar measuring device can process measured values between 20 and 630 mg/dL. The warning message "Lo" is displayed for measured values below 20 mg/dL. The warning message "Hi" is displayed for a measured value above 630 mg/dL.



#### Note

- If you suspect that the blood sugar values are incorrect, first repeat the test and, if necessary, perform a function test using the control solution. If you consistently have reason to doubt the results, contact your doctor.

#### Blood sugar

The following table shows a classification of blood sugar values according to the guidelines of the German Diabetes Association (DDG).

<b>Time of blood sugar measurement</b>	<b>Normal blood sugar values</b>	<b>Suspicious</b>	<b>Diabetes</b>
On an empty stomach	Below 100 mg/dL	100-110 mg/dL	Above 110 mg/dL
2 hours after eating	Below 140 mg/dL	140-200 mg/dL	Above 200 mg/dL

Source: Deutsche Diabetes Gesellschaft (DDG) 2008

## Evaluating critical measured values

Display	Blood sugar	Remedy
<b>Lo</b> mg/dL	Hypoglycaemia Below 20 mg/dL	Immediate medical attention re- quired.
<b>65</b> mg/dL	Low blood sugar Below 70 mg/dL	Eat a suitable snack or small meal. Follow your physician's instructions.
<b>150</b> mg/dL	High blood sugar With an empty stomach Above 100 mg/dL  2 hours after eating Above 140 mg/dL	If the value remains high 2 hours after your last meal, this may indi- cate hyperglycaemia (high blood sugar). Discuss with your doctor any possible action to take in this case.
<b>300</b> mg/dL	High blood sugar, possibly ke- tones Above 240 mg/dL	Perform a ketone test. Consult your doctor.
<b>H<sub>i</sub></b> mg/dL	Very high blood sugar Above 630 mg/dL	Repeat the measurement with a new test strip. If the display shows the same result again: Seek medi- cal advice immediately.



## 5.4 Function check with control solution

---

The control solution is used for checking the whole blood sugar measuring system. This enables you to determine whether the measuring device and the test strips are functioning correctly together and whether the test is performed correctly. It is very important that you perform these tests using the control solution to ensure accurate measuring results.



### Important

Never use the control solution of other manufacturers. You can only verify the proper function of your measuring device using the Beurer GL40 control solution.

### When is it recommended to perform a test using the control solution?

- After opening a new test strip container.
- If you suspect that the measuring device or the test strips could be defective.
- If your blood sugar measurement values do not correspond to the way you are feeling.
- If you have dropped the measuring device or it has been subjected to other mechanical strain.
- Every time you suspect that the blood sugar results may be incorrect.

It is sufficient to perform a single control solution test if the result is within the recommended range.



### Warning



Do not ingest the control solution. The control solution is only to be used for function tests and is for external use only.






### Note

- Store the control solution tightly closed at a room temperature below 30°C. Do not refrigerate.
- Store the control solution out of reach of young children.
- Do not store the control solution in the medicine cabinet. It may be mistaken for medicine to be consumed.
- The control solution can discolour your clothing. In this case, wash your clothing with soap and water.
- Carefully close the control solution bottle after each use.
- Any control solution squeezed out in excess must not be sucked back into the bottle.



### Preparations

- Do not use the control solution once the expiry date has passed. Using control solution that is out-of-date may lead to incorrect function test results. The expiry date is located on the bottle next to the hourglass symbol .
- The control solution can be kept for three months once the bottle is open. Make a note of when this time will run out (date of opening + 3 months ) on the label. The shelf

life decreases if the expiry date is reached before the end of the 3 months (see the date next to the hourglass symbol ).





- Do not use the control solution if either of these dates has expired (  /  ).
- Allow the measuring device, test strips and control solution to reach room temperature (+20°C to +25°C).

## Performing a function test with control solution

- 1 Hold the measuring device so that the display is facing you.
- 2 Insert a test strip, contacts first, into the slot on the measuring device. Ensure that the front side of the test strip is facing you (see "Test strips" page 14).
- 3 The unit switches on automatically and briefly shows the start display. As soon as the hand  and the  symbol are flashing, the unit is ready to perform measurements.



### Note

Press the  or  rocker button to switch to control mode. "EEL" appears in the display. This means that the result is not saved in the memory so that your measurement statistics are not falsified. Pressing  or  again causes "EEL" to disappear from the display again and the value will be stored as usual in the memory.

- 4 Shake the control solution well before use. Unscrew the cap and squeeze out a drop of solution. Wipe away the first drop and squeeze out an additional drop.
- 5 To prevent the remaining control solution in the bottle from becoming contaminated through the tip of the bottle in contact with the test strip, do not apply the drop directly to the test strip. Apply the drop to a clean substrate. Then put the drop on the slit for receiving blood on the test strip. The solution is sucked into the slit. Wipe the tip of the bottle using a clean, dry paper tissue.
- 6 Once the slit is filled with solution, the unit begins the measurement. The unit counts down for approx. five seconds. The result of the measurement is then shown in the display.
- 7 Check whether the result is within the specified results range for the control solution. This results range is printed on the container with the test strips.

## Results to be expected

At room temperature, the measurement results for 95% of all control solution tests should lie within the results range printed on the test strip container.

### **Warning**

The results range printed on the container with the test strips applies only to the control solution. **This is not a recommended value for your blood sugar content.**

If measurement results are outside the specified range, check the following possible causes:

Cause	Remedy
<ul style="list-style-type: none"><li>• The first drop of control solution was not disposed of.</li><li>• The tip of the bottle was not wiped clean.</li><li>• The bottle was not shaken vigorously enough.</li></ul>	Eliminate the cause and repeat the test.
The control solution or the test strip has expired or is contaminated.	Repeat the test with a new bottle of control solution or with a new test strip.
The control solution, test strip or measuring device is too warm or too cold.	Allow the control solution, test strips and measuring device to reach room temperature (+20°C to +25°C) and repeat the test.
Test strip is damaged.	Repeat the test with a new test strip.
Test strip is out of date.	Open a new container of test strips. Repeat the test.
There is a problem with the measuring device.	Please contact customer service.

### **Warning**

If you repeatedly obtain measurement results with the control solution which are outside of the specified range, **refrain from using the system to determine your blood sugar content.** Please contact customer service.

## 6 The Measured Value Memory

Your blood sugar level is automatically stored with every measurement, unless "CCL" was activated during a blood sugar measurement with control solution.

The measured value memory can store up to a maximum of 480 values. Once this number is reached, the measured value always replaces the oldest value in the memory. You can retrieve each individual blood sugar measured value. For the blood sugar levels, you can calculate and display the average value for the last 7, 14, 30 and 90 days.

### Note

- If measured values are already stored in the memory and you then reset the date, the average values are calculated according to the new time period.
- "----" indicates that the measured value memory is empty. Press the ON/OFF button to switch off the unit.

### 6.1 Displaying individual values

The individual values of the last 480 measurements are displayed. The most recent value is displayed first, and the oldest value is displayed last. The measuring device also displays the date and time of the measurement.

1 The measuring device must be switched off. Press the ▲ or ▼ rocker button [5].

2 The start display appears briefly.

"MEM" and the number of stored blood sugar tests are displayed briefly (Fig. 1). Then the stored value is displayed with the measuring unit, date, time and "MEM" (Fig. 2).





Fig. 1

Fig. 2

3 Every time you press the ▼ rocker button [5] again, the memory position number appears followed by the previous measured value. You can display a maximum of 480 previous measured values.

4 You can cancel the process at any time. To do so, press the ON/OFF button or wait until the unit switches off automatically after 1 minute.


## 6.2 Displaying average blood sugar levels

---

You can display the average blood sugar value from the last 7, 14, 30 and 90 days.

- 1 The measuring device must be switched off. Press the ▲ or ▼ rocker button [5] for 3 seconds.  
The start display appears briefly.  
The blood sugar level unit, "7 days" and the average value are displayed.  
The number of stored measured values is displayed. Press the ▲ rocker button [5] again.
- 2 Press ▲ [5] repeatedly to display the average value for 7, 14, 30 and 90 days.
- 3 You can cancel the process at any time. To do so, press the ON/OFF button or wait until the unit switches off automatically after 1 minute.

Item	Meaning
1	Average value
2	Number of days, e.g. 7



## 6.3 Deleting the measured value memory

---

- 1 The measuring device must be switched off.
- 2 Remove the batteries from the battery compartment.
- 3 Simultaneously press the ▲ button [5] and the ON/OFF button [3] and hold them down.
- 4 Reinsert the batteries while holding down the buttons; pay attention to the correct polarity.
- 5 After inserting the batteries, continue holding down the ▲ button and the on/off button [3] for another five seconds. "dEL" is displayed and the measured values are deleted.

## 6.4 Transferring measured values to a computer

---

The GL40 measuring system has a built-in computer interface [4] which you can use to transfer the measured values stored in your unit to a computer. (For the position of the connection socket, see p. 12).

You can acquire the connection cable as the accessory set 'Beurer GL40 PC Kit' in a specialist store (see "Delivery scope, replacement parts and accessories"). Enclosed with the connection cable is a CD with 30-day testing software for evaluating your measurement results, which makes it easier for you and your physician to follow your blood sugar trends.

For more information, please refer to the instructions included with the accessory set. These contain all the information you require for the data transfer.



### Note

- An effective evaluation is possible only if you have correctly set the date and time (see p. 17).
- It is not possible to take a measurement during the data transmission.
- The measurement data remain stored on the measuring device after being transferred to the computer.



### Important


Only use the original data connection cable supplied by Beurer for data transfer. Otherwise, your measuring device or PC may be damaged.

## Preparations

- Set the blood sugar measuring device near your computer.
- Install the evaluation software on your computer as outlined in the instructions for the Beurer PC Kit.

## Transferring your measurement data

- 1 The measuring device must be switched off. Insert the larger, flat USB connector of the connection cable into a USB port on your computer. Insert the mini USB connector into the computer interface port on your measuring device.
- 2 "-> PC" appears in the display of the measuring device. The measuring device is now ready for data transfer.
- 3 Follow the transferring and evaluating information in the software and in the instructions for the Beurer PC Kit.



## 7 Storage and maintenance of the unit

### Storing your unit

Store the measuring device after each use in the case provided.



#### Note

- Keep these instructions for use.
- If you will not use the unit for a prolonged period, remove the batteries.

### Maintenance

The surface of the unit can be cleaned with a damp cloth (moistened with water or a mild cleaning solution). Dry the unit using a lint-free cloth.



#### Note

The measuring device is made from precision components. The accuracy of the measured values and the service life of the unit depend on careful handling:

- Protect the unit against impact and do not drop it.
- Protect the unit against damaging influences such as moisture, dirt, dust, blood, control solution or water, extremes of temperature and direct sunlight.
- Do not store close to strong electromagnetic fields, wireless systems or mobile telephones.



## 8 Troubleshooting



### Messages on the display for batteries and blood sugar measurement

No.	Cause	Action
LP	Flat batteries	Replace all batteries.
Ht	Temperature of the measurement room, measuring device or test strip above the allowable limit.	Repeat test with new test strip as soon as the measurement room, measuring device and test strip have reached room temperature (+20°C to +25°C).
Lt	Temperature of the measurement room, measuring device or test strip below the allowable limit.	Repeat test with new test strip as soon as the measurement room, measuring device and test strip have reached room temperature (+20°C to +25°C).
Err	Used or contaminated test strip was inserted.	<ul style="list-style-type: none"><li>• Insert test strip which has not been used or spoiled.</li><li>• Repeat blood sugar measurement.</li></ul>
001	System error	<ul style="list-style-type: none"><li>• Remove the batteries and reinsert them. Contact our customer service department if the problem persists.</li></ul>
005	Memory error	<ul style="list-style-type: none"><li>• Remove the batteries and reinsert them. Contact our customer service department if the problem persists.</li></ul>

**Problem: Unit does not turn on.**

Cause	Action
Flat batteries	Replace the batteries.
Incorrectly inserted or missing battery.	Verify that the batteries are inserted correctly (see "Inserting and replacing the batteries", page 16).
Test strip is inserted with the wrong side or not completely.	Insert the test strip, contacts first, into the slot on the unit. Ensure that the front side of the test strip is facing you (see "Test strips", page 14).
Unit is defective	Contact customer service.

**Problem: After inserting the test strip into the unit and applying the blood, the test does not start.**

Cause	Action
Insufficient quantity of blood	Repeat the test using a new test strip and larger blood drop.
Test strip is defective	Repeat the test using a new test strip.
Blood was applied when the unit was switched off.	Repeat the test and apply the blood only when  and  are flashing.
The unit's basic settings were changed and the change was not completed (see "Basic settings", page 17).	Press the "On/Off" button until "OFF" is displayed. Repeat test.
Unit is defective	Contact customer service.

## 9 Technical specifications

Dimensions (WHD)	47 x 85 x 14 mm
Weight	43 g (including batteries)
Electrical power	2 x 3V CR2032 button cell batteries
Battery life	More than 1000 measurements
Measured value memory	480 measured values with date/time Data retention when replacing the batteries
Average values	Blood sugar: for 7, 14, 30, 90 days
Automatic switch-off	2 minutes after the last operation
Storage/ transport temperature	Temperature: +2°C – +30°C Relative humidity: < 90%
Operating ranges	Temperature: +10°C – +30°C Relative humidity: < 90% non-condensing
Glucose measurement range	Glucose: 20 – 630 mg/dL
Blood sample	Capillary whole blood
Blood quantity	0.6 microlitres
Measuring time for blood sugar	Approx. 5 seconds
Calibration	Plasma
Test method	Amperometric biosensor
Application	Suitable for self testing
System function test	Each time the unit is switched on

## **EMC**

This unit complies with European standard EN 61326 and is subject to specific safety precautions in terms of electromagnetic compatibility. Note that portable and mobile RF communication equipment can affect this unit. Further details can be requested by contacting the customer service address provided below.

### **Note regarding function of test strips**

The test strips enable a quantitative measurement of the glucose content of capillary whole blood. When the slit for collecting blood comes into contact with a drop of blood, it fills automatically through simple capillary action. The blood is sucked into the absorbing slit of the test strip and the measuring device measures the sugar level of the blood.

The test is based on the measurement of an electric current caused by a chemical reaction between the glucose and the test strip reagent. The measuring device analyses this current. The flow of current depends on the glucose content of the blood sample. The results are shown in the display of the measuring device. Only a small quantity of blood is required (0.6 microlitres) and the measurement takes approx. five seconds. The test strips can detect blood sugar levels in the range from 20 to 630 mg/dL.

### **Chemical components of the test strip sensor**

- Glucose oxidase  $\geq 0.6$  IU
- Electron shuttle  $\geq 0.03$  mg
- Non-reactive components  $\geq 0.06$  mg

### **Note regarding function of the control solution**

The control solution contains a fixed concentration of glucose that reacts with the test strip. A test with the control solution is performed in the same way as a blood test, but the control solution is used instead of blood. The measurement result from the control solution must lie within the results range. This results range is printed on every test strip container.

### **Chemical composition of the control solution**

The control solution is a red-coloured solution with a D-glucose content of less than 0.2%.

Ingredient	Percentage concentration
D-glucose	0.12 %
Non-reactive components	99.88%

### **Standards**

The Beurer GL40 measuring system for household use corresponds to the following European directives and standards:

IVD (98/79/EC), EN 61010-1, EN 61010-2-101, EN 13640, EN ISO 15197, MDD (93/42/EC).

In accordance with the MPBetreibV (German Medical Device Operator Ordinance), regular metrology checks must be performed if the unit is to be used for commercial or business purposes. For private use, we also recommend having the manufacturer perform a metrology check every 2 years.

### **Comparison between measured values and laboratory values**

Performance characteristics: Accuracy and precision

Whole blood sugar test results were compared with the YSI 2300 laboratory instrument. For a concentration < 75 mg/dL, ≥ 95% lay within +/- 15 mg/dL, while for a sugar concentration ≥ 75 mg/dL, ≥ 95% lay within 20% of the reference values. The CV (variation coefficient) is < 5 %. The blood sugar measuring device is therefore comparable with a laboratory system.

## 10 Guarantee and Customer Service

### Guarantee

This product comes with a 3-year guarantee for material and manufacturing faults. The guarantee does not apply:

- in the case of damage caused by improper use
- to wearing parts
- to deficiencies of which the customer was aware at the time of purchase
- to personal negligence on the part of the customer
- in the case of third-party intervention

This guarantee does not affect your statutory rights. In order to make a claim within the warranty period, the customer is required to provide proof of purchase. Claims must be made within a period of 3 years from the date of purchase to BEURER GmbH, Söflinger Str. 218, 89077 Ulm, Germany. In the case of claims against the guarantee, the customer has the right to have the product repaired by us or in a workshop authorised by us. Further rights (of the guarantee) remain unaffected.

### Customer service address


If you have any questions, please contact customer service:

Ⓒ Lifestyle MI Ltd  
P.O. Box 584  
WN1 9 EX WIGAN  
Phone +44 870 879 08 12  
E-Mail: [customerservice@lifestylemi.com](mailto:customerservice@lifestylemi.com)

Ⓘ Brandlinx Direct Limited  
Hainault House  
Baldonnell Duriken Park  
Dublin 22  
Phone +353 1 412 3606  
E-Mail: [sales@brandlinx.ie](mailto:sales@brandlinx.ie)

**OUR COMMITMENT TO YOU:** We aim to satisfy our customers by providing high-quality healthcare products and the best customer service. If you are not completely satisfied with this product, please contact customer service.





GL40\_mg-dL\_1208\_GB Subject to error and change