BELL & HARP

AUDIOMETERS

QUICK START GUIDE

EN – Quick Start Guide

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QUICK START GUIDE: ENGLISH

INTENDED USE

Bell & Harp audiometers are designed for use by audiology and hearing healthcare specialists or trained technicians in conducting diagnostic hearing evaluations and assisting in the diagnosis of possible otologic disorders. Audiometric tests must be conducted in particularly quiet environment, and preferably in a soundproof booth. If these conditions cannot be provided, the patient should be given sound-isolating headphones to wear, so that ambient and background noise can be blocked out as far as possible.

PRECAUTIONS

Users are entirely responsible for any malfunction caused by improper use, or by maintenance or repair operations performed by any party other than Inventis srl or an authorised Service Centre.

Installation

If the audiometer is connected to other devices, obtaining an electromedical system, this system must be complying with IEC standard 60601-1.

Avoid installing and using Bell & Harp audiometer close to any sources of strong electromagnetic field: these could interfere with the operation of the device.

Calibration

The calibration is valid only for the transducers supplied with the device. If a transducer is replaced, the audiometer must be recalibrated.

Use

Make attention to set a suitable intensity of stimulus before presented it to the patient.


The eartips of ER-3 and ER-5 insert earphones and the masking earphones are disposable: do not use the same for different patients. Dispose of eartips after use.


Disinfect the cushions of headphones between one patient and the next.


REAR PANEL CONNECTIONS

**Rear Panel of Harp audiometer.
Bell audiometer does not have some of these connectors.**

<i>Connector</i>		<i>Description</i>
1		Power switch (0: Off, 1: on)
2		Power supply connection
3	USB	USB cable for the connection to the PC
4	PATIENT RESP. 	Patient response switch
5	TALK OVER	Microphone for operator to patient communication
6	MONITOR	Monitor headphones
7	EXT	External input 1 for speech audiometry
8		External input 2 for speech audiometry
9	MIC	Operator microphone for live speech tests
10	TALK BACK	Patient microphone
11	SPKR	Free field speaker: Left L
12		Free field speaker: Right R
13	PHONES 	AC headphones - Right R
14		AC headphones - Left L
15	INSERT 	Insert Masking earphone
16	BONE VIBR. 	Bone vibrator
17	INSERT	AC Insert earphones - Right R
18		AC Insert earphones - Left L

THE KEYBOARD

<i>Switch</i>	<i>Function</i>
1	Starts or pauses the playback of words list.
2	Terminates playback of the current list.
3	If audio material is not playing: brings up the list preceding the current list. If audio material is playing: allows to jump back to the previous word, if the list is indexed, or by 5 seconds if not indexed.
4	If audio material is not playing: brings up the list that follows the current list. If audio material is playing: allows to jump forward to the next word, if the list is indexed, or by 5 seconds if not indexed.
5	BACK Allows to go back to the previous window.
6	NEW PATIENT Cancels all examinations conducted in the current session.
7	SAVE EXAM Saves the current test to the patient memory of the instrument, associating it with the patient indicated in the main window.
8	ERASE DATA Deletes the current data.
9	TO PC PRINT Sends the current examination to a computer. Starts the printout of the current test (if the audiometer is equipped with the thermal printer).
10	HELP Opens a context-sensitive help window.
11	dB Selects the incremental step for the attenuators between 1 or 5 dB.

12	PULSED	Allows the selection of pulsed mode for the stimulus signal.
13	TRACKING	If enabled, the difference between the intensity of the stimulus and masking signals is maintained constant, changing the intensity of the stimulus signal.
14	LOCK	If enabled, the masking tone will be switched on and off with the stimulus tone (The interrupter and the Normally ON button of Channel 2 are deactivated).
15	TALK	Shows a volume control window for communications between patient and operator.
16	GAIN	(Only for Speech audiometry) Shows a gain control window for the speech material.
17	MONITOR	Shows a volume control window for monitor signals used by the operator.
18	RIGHT	Selects the stimulation side for the Channel 1. For stimulation on both sides, press the two keys simultaneously.
19	LEFT	
20	OUTPUT	Selects the transducer for the Channel 1.
21	INPUT	Selects the input signal for the Channel 1.
22	INPUT	Selects the input signal for the Channel 2.
23	OUTPUT	Selects the transducer for the Channel 2.
24	RIGHT	Selects the stimulation side for the Channel 2. For stimulation on both sides, press the two keys simultaneously.
25	LEFT	
26	NORM. ON	If enabled, the output of Channel 1 is permanently active. The channel is switched off by touching the interrupter (31).
27	TALK OVER	Enables, by holding the button, the operator-to-patient communication.
28	TALK BACK	Enables the patient-to-operator communication.
29	NORM. ON	If enabled, the output of Channel 2 is permanently active. The channel is switched off by touching the interrupter (32).
30	STORE	In pure-tone audiometry: saves the threshold at each frequency. In all other tests: saves the values measured.
31	INTERRUPTER	Switches on the Channel 1 (31) or Channel 2 (32).
32		

33	FREQUENCY	Changes the frequency value to the previous and next one. In the case of Speech audiometry: allows to count the words recognized (CORRECT) and not recognized (INCORRECT) by the patient.
34		
35	Knob	Allows to attenuate the output or scrolls a list of settings.
36		
	RESPONSE	Lights up whenever the patient presses the response button.

OPERATING INSTRUCTION

* To select those functions marked with the asterisk, press the corresponding key, placed below the display.

Selecting an exam

Scroll the list of exams with the left knob (35) and open the specific Menu by pressing *OK** or (31).

Communicating with patient

Press and hold (27). Communicate through the microphone built into the chassis or through an external microphone plugged into connector (5). Adjust the volume by (35).

Pure-tone audiometry

- Select the input signal for Channel 1 (stimulus): (21);
- Select the input signal for Channel 2 (masking): (22);
- Select the transducer for Ch.1 (20) and for Ch.2 (23);
- Select the mode of presenting the stimulus (12) between continuous or pulsed and, in case, change the pulse frequency by pressing the appropriate key below the display;
- Select the stimulation side: R (18), L (19) or Bilateral (18+19); the masking signal is automatically set to contralateral;
- Select the frequency (33), (34) and the intensity (35) of the stimulus; to select an intensity above 100 dB HL, press the *HIGHER dB** button;
- In case you need the masking, select the intensity (36) and enable it (29);
- Deliver the stimulus (31);
- When the threshold is identified, store it (30) or press *NO RESP.** if the patient cannot hear the tone at its maximum intensity.

Auto Threshold

- Select the input signal for Channel 1 (stimulus): (21);
 - Select the input signal for Channel 2 (masking): (22);
 - Select the transducer for Ch.1(20) and for Ch.2 (23);
-

- Select the side of stimulation (18), (19) and of masking (24), (25);
- Set the frequency to test: open the relevant window by *FREQ.SELECTION** and select the frequency value using (35), (36);
- Start the test, by pressing *START EXAM**.

Speech Audiometry

- Select the input signal for Channel 1: (21); the INT input refers to the internal flash memory;
- Select the masking signal: (22);
- Adjust the inputs gains: play the calibration tone (if the input is EXT or INT), press (16) and regulate the gains through (35), (36);
- Select the transducers for Ch.1(20) and for Ch.2 (23);
- Select the stimulation side: R (18), L (19), or Bilateral (18+19); the masking signal is automatically set to contralateral;
- Select the intensity of stimulus (35) and masking signal (36);
- Play the words list (in case of INT input, use (1), (2), (3) and (4)) or say the words (in case of MIC input);
- Update the score, by pressing (33) or (34), depending on the Speech score mode and the correct/incorrect patient's answer;
- To save the score, press (30).

For uploading a words list onto the internal flash memory, use the *Audio Tracks Indexing Tool (ATIT)* software, supplied with the audiometer. With this software you can also *index* the lists, so that you can see on the device display the word being pronounced. See the Bell & Harp User Manual for a detailed description of ATIT.

Supra Threshold tests, Stenger and MHA

ABLB (Fowler)

- Set the side - (18), (19) or (24), (25) - to associate to Channel 1(*Ref*: the good ear) and to Channel 2(*Var*: the poor ear);
 - Select the frequency to test (33), (34) and adjust the signal intensity for each channel: (35), (36);
 - Select the frequency for alternating the tone between two ears (*PULSE RATE**);
 - Present the signals (31);
 - Save the point in which the patient perceives the signals as equals: (30).
- Once four tests are saved, navigate among these using (33), (34).

SISI (Jerger)

- Select the masking signal on Channel 2: (22);
 - Select the transducers for Ch.1 (20) and for Ch.2 (23);
 - Select the side to test: (18), (19);
-

- Set the value of tone frequency (33), (34) and adjust the starting intensity value of tone (35) and masking (36);
- Set the amplitude of increments, by pressing (11) and the rate for presenting these increments (*PULSE RATE*);
- Start the test, by pressing *START EXAM**.

DLI (Lüscher)

As the SISI test, with the only difference the operator can modify the amplitude of increments through *dB INCREM.-* and *dB INCREM.+* and have to save the level of the minimum increment of intensity identified by the patient, by pressing (30).

Stenger

- Select the input signal: (21); the input can be either a tone or a speech signal;
- Adjust the intensity value of stimulus for right (35) and left (36) ear and, in case of a tone, set the frequency to test (33), (34);
- Present the signal: (31).

No results are saved for this type of test.

Tone decay

- Set the duration of the exam, by pressing *DURATION**;
- Select the input signal for Channel 2 (masking): (22);
- Select the transducer for Ch.1 (20), and for Ch.2 (23);
- Select the stimulation side: R (18), L (19);
- Select the frequency (33), (34) and the intensity (35) of the stimulus;
- In case of masking, select the intensity (36) and enable it (29);
- Start the test, by pressing *START EXAM**.

Master Hearing Aid

- Set the side for Channel 1 (18), (19) or (18+19) and Channel 2 (24), (25) or (24+25);
- Select the input signal (21);
- Set the filter for the Channel 1 (*FILTER CH.1**) and for the Channel 2 (*FILTER CH.2**);
- Adjust the intensity value of stimulus: (35), (36);
- Turn on the channel(s): (26), (29) and start presenting the speech material.

No results are saved for this type of test.

TECHNICAL SPECIFICATIONS

These technical specifications cover the main aspects of the device: more specific details are provided in the User Manual.

Classification

Bell Basic. Pure-tone audiometer: Type 4 (IEC 60645-1 / ANSI S3.6)

Bell Plus. Pure-tone audiometer: Type 3 (IEC 60645-1 / ANSI S3.6)

Harp Basic and Plus

Pure-tone audiometer: Type 2 (IEC 60645-1 / ANSI S3.6)

Speech audiometer: Type A or A-E (IEC 60645-1 / ANSI S3.6)

Channels

Two independent channels.

Available signals

Stimulus: pure tone, warble tone.

Masking: NBN, WN (not Bell Basic), SN (only Harp).

Speech audiometry inputs: EXT1, EXT2, MIC, INT (flash memory, only Harp Plus).

Signals specifications

Attenuator step: 5 dB (Bell); 1 and 5 dB (Harp).

Presentation: Continuous, Pulsed (0.5, 1, 2 Hz).

Warble: 5 Hz sin wave modulating signal

Available outputs and transducers

AC: TDH-39 or DD45 headphones, ER-3 insert earphones

BC: B-71 bone vibrator

Free field (only Harp)

Insert masking earphone: IME-100

Pure tone: available frequencies and maximum levels (dB HL)

Freq. (Hz)	AC TDH39 DD45	AC ER-3	AC ER-5	BC	FF(*)
125	80	90	90	-	75
250	100	105	100	45	85
500	110	110	110	65	95
750	115	115	120	70	95
1.000	120	120	120	75	95
1.500	120	120	120	80	95
2.000	120	120	115	80	95
3.000	120	120	115	75	95
4.000	120	110	110	75	95
6.000	115	95	100	55	90
8.000	100	90	90	50	85

(*) The values refer to “normal” range; add 10 dB to each value in case of “extended range” option selected.

Speech audiometry: maximum levels (dB HL)

	AC: <i>TDH39</i> <i>DD45</i>	AC: <i>ER-3</i> <i>ER-5</i>	BC	FF
Type A	100	100	60	Normal: 75 Extended: 85
Type A-E	80	100	60	Normal: 75 Extended: 85

Patient – operator communication and monitoring

Talk over: via built-in or external microphone.

Talk back: patient microphone input.

Available tests

AVAILABLE TESTS				
	<i>Bell</i> <i>Basic</i>	<i>Bell</i> <i>Plus</i>	<i>Harp</i> <i>Basic</i>	<i>Harp</i> <i>Plus</i>
Pure tone audiometry	●	●	●	●
Autothreshold	●	●	●	●
Speech audiometry	-	-	●	●
ABLB	-	-	●	●
SISI	-	-	●	●
DLI	-	-	●	●
Tone decay	-	-	●	●
Master hearing aid	-	-	●	●
Stenger	-	-	●	●

Internal flash memory (only Harp Plus)

Used to store the speech material (*wav* format); *Capacity:* 2 or 4 GB.

Speech material upload: through ATIT software (incl.)

Printer (only Harp)

Optional integrated thermal printer; *Paper size:* 114 mm.

Calibration

Validity: 12 months.

TDH 39, DD45: ISO 389-1; *ER-3* and *ER-5:* ISO 389-2; *B71:* ISO 389-3;

FF: ISO 389-7.

Computer interface

Connection: USB (driverless); *Compatible software:* Inventis Maestro.

Display

Type: Graphical colour TFT LCD; *Size:* diagonal 4.3", 95 mm x 54 mm;
Resolution: 480 x 272

Power supply

Consumption: 9 Watts; *Power supply:* 6V, 2A cont., through an external medical grade 100-240 Vac 50/60 Hz power supply.

Mechanics*Without integrated printer*

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