ARGON MEDICAL DEVICES, INC.
1445 Flat Creek Road, Athens, Texas 75751 USA
Tel: +1 (903) 675-9321;
Tel: +1 (800) 927-4669
www.argonmedical.com

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands
+31 70 345 8570

Argon Medical Devices UK Ltd
Eastgate Business Centre
Eastern Avenue
Burton-on-Trent
Staffordshire
DE13 0AT

IFU9011/0321A

<https://www.argonmedical.com/resources/product-information>

English	The symbols glossary is located electronically at www.argonmedical.com/symbols
Bulgarian	Речникът за символите в електронен вариант се намира на www.argonmedical.com/symbols
Czech	Slovníček symbolů je elektronicky umístěn na www.argonmedical.com/symbols
Danish	Symbolforklaringen findes elektronisk på www.argonmedical.com/symbols
German	Das Glossar der Symbole kann unter www.argonmedical.com/symbols abgerufen werden
Greek	Το γλωσσάρι των συμβόλων παρέχεται ηλεκτρονικά στη διεύθυνση www.argonmedical.com/symbols
Spanish	El glosario de símbolos puede consultarse electrónicamente en www.argonmedical.com/symbols
Estonian	Sümbolite seletuste kogu on elektrooniliselt saadaval aadressil www.argonmedical.com/symbols
Finnish	Elektroninen symbolisanasto on osoitteessa www.argonmedical.com/symbols .
French	Le glossaire des symboles est disponible au format électronique à l'adresse www.argonmedical.com/symbols
Croatian	Pojmovnik simbola nalazi se u elektroničkom obliku na www.argonmedical.com/symbols
Hungarian	A szimbólumok összefoglalása elektronikusan itt található: www.argonmedical.com/symbols
Italian	Il glossario dei simboli è disponibile in formato elettronico all'indirizzo www.argonmedical.com/symbols
Lithuanian	Simbolių žodyno elektroninė forma pateikta www.argonmedical.com/symbols
Latvian	Simbolu glosārijs ir atrodams elektroniskā veidā vietnē www.argonmedical.com/symbols
Dutch	Uitleg over de symbolen is elektronisch beschikbaar op www.argonmedical.com/symbols
Norwegian	Symbolforklaringen befinner seg elektronisk på www.argonmedical.com/symbols
Polish	Glosariusz symboli w postaci elektronicznej dostępny jest pod adresem www.argonmedical.com/symbols
Portuguese	O glossário de símbolos está localizado eletronicamente no endereço www.argonmedical.com/symbols
Romanian	Glosarul simbolurilor este stocat electronic la www.argonmedical.com/symbols
Slovak	Glosár symbolov je v elektronickej podobe k dispozícii na stránkach www.argonmedical.com/symbols
Slovenian	Slovar simbolov se nahaja v elektronski obliki na www.argonmedical.com/symbols .
Swedish	Symbolforklaringen hittas elektroniskt på www.argonmedical.com/symbols
Turkish	Semboller sözlüğü elektronik ortamda www.argonmedical.com/symbols adresinde yer almaktadır

EN - ENGLISH

Techna-Cut™ Biopsy Needle

Intended Use/Purpose

Techna-Cut™ Biopsy Needles are intended to be used with a range of biopsy application of soft tissue.

Device Description

The Techna-Cut™ Biopsy Needle is a cutting biopsy needle with a trocar point stylet and sharp cannula. A locking syringe is included with the needle. The product are available in 16ga, 18ga, 19ga, 21ga and 23ga in lengths of 6cm, 7cm, 12cm, 10cm and 15cm.

Indications for Use:

The Techna-Cut™ Biopsy Needle is intended for use in obtaining biopsy material from tissues such as the liver, kidney, prostate, breast, thyroid, pancreas, spleen, and lung.

Warnings

- The device is designed, tested and manufactured for single use only.
- Reuse or reprocessing has not been evaluated and may lead to its failure and subsequent patient illness, infection or other injury.
- Do not reuse, reprocess or re-sterilize the device.
- Inspect the package integrity before use.
- Do not use if package is open or damaged and if the expiry date has been exceeded.
- Do not attempt to use the device on hard tissue.
- Do not use excessive force, it may result in damage/breakage of the device.
- Do not bend the needle excessively, it may result in damage/breakage of the device.
- Verify the integrity of the needle before inserting. If needle is damaged, replace the needle.

Precautions:

- These devices should be used by a physician familiar with the possible side effects, typical findings, limitations, indications and contraindications of needle biopsy.

Possible Complications:

Bleeding
Pain
Swelling at site
Infection

NOTE: These instructions are NOT meant to define or suggest any medical or surgical technique. The individual practitioner is responsible for the proper procedure and techniques to be used with this device.

BIOPSY PREPARATION/PROCEDURE:

- Using sterile technique, prepare the biopsy puncture site.
- Insert the needle through the skin with the stylet in place, and advance the needle, typically under imaging guidance, directly into the lesion.
- Remove the stylet and attach a syringe.
- While applying suction, move the needle up and down rapidly several times within the lesion.
- When the aspiration procedure is complete, maintain gentle suction and withdraw the needle.
- The sample (aspirate or tissue fragments) are now contained within the needle and ready to be deposited on a slide or in formalin for cytological analysis.
- Remove the needle.
- Apply pressure to puncture site to control possible bleeding.

Disposal

After use, this product may be potential biohazard. Handle in a manner, which will prevent accidental sharps puncture. Dispose in accordance with hospital policies and procedures for waste and biohazard materials.

Storage

Store in ambient temperature.

NOTE: In the event a serious incident related to this device occurs, the event should be reported to Argon Medical at quality.regulatory@argonmedical.com as well as to the competent health authority where the user/patient resides.

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