

Une chaleur apaisante

La chaleur pénètre en profondeur grâce au diffuseur de chaleur infrarouge 300 W. Elle favorise la circulation sanguine et contribue au processus de guérison et à la détente des muscles et des tensions musculaires.

Soothing warmth

The deeply penetrating warmth of the 300 W infrared radiator promotes blood circulation and can therefore help to support the healing process and relax muscles and muscle tension.



germanengineering

IL 50 Diffuseur de chaleur infrarouge

Vitre en verre céramique

Design exclusif

Angle d'inclinaison réglable en continu

300 watt / 220-240 V 50-60 Hz

Champ d'irradiation 30 x 40 cm

Taille : 27 x 28,5 x 19,5 cm

Poids : 1,6 kg

2 m câble

Enroulement du câble

Changement des tubes possible

Minuteur réglable de 1 à 15 min ;
arrêt automatique

Avec ventilateur

Sécurité anti-surchauffe

Produit médical

Garantie 3 ans

Colisage : 2

EAN-Nr.: 42 11125 61900 8

Art.-Nr.: 619.00

Carton de transport : -

IL 50 Infrared heat radiators

Ceramic glass

Exclusive design

Infinitely adjustable shade

300 watt / 220-240 V 50-60 Hz

Area of radiation 30 x 40 cm

Measurement : 27 x 28.5 x 19.5 cm

Weight: 1.6 kg

2 m cord

Cord winder

Exchange of tubes possible

Timer adjustable 1-15 min;
with automatic switch off

With fan

Overheating protection

Medical product

3 year guarantee

Sales unit: 2

EAN no.: 42 11125 61900 8

Item no.: 619.00

Shipping carton: -