

PRODUCT LIST 2020



The innovative way of Cryosurgery



Phone ► +46 31-400 500
Fax ► +46 300-56 94 99

E-mail ► order@cryoiq.com
Website ► www.cryoiq.com



How to order

If possible please send your order via email. For customers within the EC community it is important that the VAT number is included in the order.

Order by e-mail: order@cryoiq.com
Order by Phone: +46 31-400 500
Order by Fax: +46 300-56 94 99

Product information

To find out more about our products, please visit our website www.cryoiq.com.

As a CryoIQ Dealer, we will also provide you with login details to our picture database. Here you will find pictures, artworks, training material, pricelist, etc. all that you need for your marketing of CryoIQ products. Please contact our sales department for your login details.

Product number configurations

Learn about product number configuration to make it easier to understand.

CIQ - AA - BB / CC - VET

CIQ CryoIQ Prefix
A Device model
- Divider between device and accessories
B Main Tip / Probe
/ Divider between extra Tip / Probe
C Extra Tip / Probe
VET Veterinary version*

* Vet version has a different UserGuide and labeling

Warranty

There is a 1+1 year warranty against technical defects. The warranty is strictly limited to the exchange of defective parts. Dropped devices or broken tips are not under warranty. Other claims for damages such as loss of working hours, incorrect treatment and their consequences, non-executed after treatment and their consequences as well as the non-observance of the safety instructions are excluded from all warranty and liability.

Warranty claims:

Since all our products are certified according to CE and registered with FDA, all warranty claims are to be registered according to specific regulations. All Warranty claims must follow the "Warrant Claim" form that can be downloaded from <http://www.cryoiq.com/warranty>. Please notice that the end-user invoice has to be attached to the Warranty Claim.

Calculating your net price

Use the attached discount table to calculate your NET prices. Note that the products are grouped into 4 discount groups A, B, C, D. You will find the corresponding groups and discounts in the discount table.

Production



CryoIQ global
Apelrödsvägen 1
439 32 Onsala
Sweden

Phone: +46 31-400 500
Fax: +46 300-56 94 99
E-mail: info@cryoiq.com
Website: www.cryoiq.com



Devices and options

CryoIQ devices are divided in 2 groups: DERM plus, PRO depending on functionality and how they could be combined with the various accessories.

DERM plus device options

Order no.	Product description	Pack. size	Discount group
STANDARD - Devices with fixed tip and gas cartridge with valve			
CIQ-DP-L	New! CryoIQ DERM plus Liquid device with 1 device for spray freezing, Ø1-6 mm 1 gas cartridge N ₂ O 25 gram 1 freeze limiter 1 case for storage 1 user guide A6 size	1	A
CIQ-DP-C1	New! CryoIQ DERM plus Contact 1 device with 1 device for contact freezing Ø1 mm 1 gas cartridge N ₂ O 25 gram 1 case for storage 1 user guide A6 size	1	A
CIQ-DP-C3	New! CryoIQ DERM plus Contact 3 device with 1 device for contact freezing Ø3 mm 1 gas cartridge N ₂ O 25 gram 1 case for storage 1 user guide A6 size	1	A
CIQ-DP-C5	New! CryoIQ DERM plus Contact 5 device with 1 device for contact freezing Ø5 mm 1 gas cartridge N ₂ O 25 gram 1 case for storage 1 user guide A6 size	1	A
CIQ-DP-L/C3/C5	New! CryoIQ DERM plus Multi pack with 1 device for spray freezing, Ø1-6 mm 1 device for contact freezing Ø3 mm 1 device for contact freezing Ø5 mm 1 freeze limiter 1 gas cartridge N ₂ O 25 gram 1 case for storage 1 user guide A6 size	1	A

DERM - Spray freezing



DPST

DERM - Contact Freezing



DPC1



DPC3



DPC5



Gas cartridge 16/25 gram



16/25 gram gas cartridge for CryoIQ for devices DERM, DERM plus and PRO.


Accessories for DERM plus devices

Order no.	Product description	Pack. size	Picture	Discount group
Replacement devices				
CIQ-DAP-DST	New! Device for spray freezing, Ø1-6 mm	1	DPST	C
CIQ-DAP-T/C1	New! Device for contact freezing Ø1 mm	1	DPC1	C
CIQ-DAP-T/C3	New! Device for contact freezing Ø3 mm	1	DPC3	C
CIQ-DAP-T/C5	New! Device for contact freezing Ø5 mm	1	DPC5	C
CIQ-DAP-T/C1-5	New! Device contact SET, Ø1, Ø3, Ø5 mm	1		C
Accessories				
CIQ-DA-C	Case for storage DERM & gas cartridges	1		C
CIQ-A-LED	LED light, for 16 gram gas cartridge	1	LED	A

Packed in white case



User's Guide A6 size



User's Guide 12 pages, with easy guidelines how to perform treatments and information on clinical studies.



Phone ► +46 31-400 500
Fax ► +46 300-56 94 99

E-mail ► order@cryoiq.com
Website ► www.cryoiq.com



PRO device options

Order no.	Product description	Pack. size	Discount group
PRO - Devices with exchangeable tip and gas cartridge with valve			
CIQ-P-L	CryoIQ PRO liquid device with 1 standard tip 60 µm, L=10 mm, spray freezing Ø1-6 mm 1 PRO body with lever 1 gas cartridge N ₂ O 25 gram New! 1 freeze limiter 1 Case for storage 1 user guide A6 size	1	A
	BEST SELLER		
CIQ-P-L/D2	CryoIQ PRO Liquid Derm 2 device with 1 standard tip 60 µm, L=10 mm, spray freezing Ø1-6 mm 1 dermatology tip D2, L=22 mm, spray freezing Ø7-18 mm 1 PRO body with lever 1 freeze limiter 1 gas cartridge N ₂ O 25 gram New! 1 case for storage 1 user guide A6 size	1	A
	BEST SELLER		
CIQ-P-C5	CryoIQ PRO Contact 5 device with 1 contact tip for contact freezing Ø5 mm 1 PRO body with lever 1 gas cartridge N ₂ O 25 gram New! 1 case for storage 1 user guide A6 size	1	A
CIQ-P-L/C5	CryoIQ PRO Mix device with 1 standard tip 60 µm, L=10 mm, spray freezing Ø1-6 mm 1 contact tip for contact freezing Ø5 mm 1 PRO body with lever 1 freeze limiter 1 gas cartridge N ₂ O 25 gram New! 1 case for storage 1 user guide A6 size	1	A
CIQ-P-L/D2/C5	CryoIQ PRO DermMix device with 1 standard tip 60 µm, L=10 mm, spray freezing Ø1-6 mm 1 dermatology tip D2, L=22 mm, spray freezing Ø7-18 mm 1 contact tip for contact freezing Ø5 mm 1 PRO body with lever 1 freeze limiter 1 gas cartridge N ₂ O 25 gram New! 1 case for storage 1 user guide A6 size	1	A
CIQ-P-L/GYN	New! CryoIQ PRO Gynecology device with 1 dermatology tip D2, L=22 mm, spray freezing Ø7-18 mm 1 gynecology tip, 5°, L=130 mm, spray freezing Ø9-22 mm 1 PRO body with lever 1 gas cartridge N ₂ O 25 gram 1 case for storage 1 user guide A6 size	1	D

PRO - Liquid Spray Freezing




PRO - Contact Freezing




User's Guide A6 size

User's Guide 12 pages, with easy guidelines how to perform treatments and information on clinical studies.



Medical version



Veterinary version



CIQ-P-L
CryoIQ PRO liquid with
1 cartridge 25 gr



CIQ-P-L/D2/C5
CryoIQ PRO liquid with Tip
D2 and C5, 1 cartridge 25 gr



CIQ-P-L/GYN
CryoIQ PRO liquid with Tip D2
and gyn. tip 5°, 1 cartridge 25 gr





Accessories for PRO devices

Order no.	Product description	Pack. size	Picture	Discount group
PRO - Accessories				
General				
CIQ-PA-DST60	Standard tip 60 µm, L=10 mm, spray freezing Ø1-6 mm	1	PDST	C
CIQ-PA-DST80	Standard tip 80 µm, L=10 mm, spray freezing Ø7-18 mm	1	PDST	C
CIQ-PA-BY	CryolQ PRO body with lever	1	Body	C
CIQ-PA-KL	Case for storage PRO & gas cartridge	1		C
CIQ-PA-CAP	Protective Cap for PRO glass tip	1	PCAP	C
CIQ-PA-LMT	Freeze limiter	1	LMT	C
Dermatology (spray freezing)				
CIQ-PA-D1	Dermatology tip 60 µm, L=22 mm, spray freezing Ø1-6 mm	1	D1	C
CIQ-PA-D2	Dermatology tip 80 µm, L=22 mm, spray freezing Ø7-18 mm	1	D2	C
CIQ-PA-D3	Dermatology tip 100 µm, L=22 mm, spray freezing Ø9-22 mm	1	D3	C
CIQ-PA-D1-3	Dermatology set PRO with 3 tips (D1, D2, D3)	1		C
Contact freezing				
CIQ-PA-C1	Contact tip for contact freezing Ø1 mm	1	C1	C
CIQ-PA-C3	Contact tip for contact freezing Ø3 mm	1	C3	C
CIQ-PA-C5	Contact tip for contact freezing Ø5 mm	1	C5	C
CIQ-PA-C7	Contact tip for contact freezing Ø7 mm	1	C7	C
CIQ-PA-C1-5	Contact set PRO with 3 tips (C1, C3, C5)	1		C
Special tip				
CIQ-PA-A45	Angled tip 45°, L=8 cm, spray freezing Ø7-18 mm	1	A45	C
CIQ-PA-L13	Long tip 100 µm, L=130 mm, spray freezing Ø9-22 mm	1	L13	C
CIQ-PA-L13/5	New! Gynecology tip 100 µm, 5°, L=130 mm, spray freezing Ø9-22 mm	1	L13/5	C
Other				
CIQ-A-LED	LED light, for 16 gram gas cartridge	1	LED	A
CIQ-SCALE-A	New! Portabel Scale, max 500 gr, 0.1 gram	1	SC	A



LED light with a flexible holder that is EASY to Snap-On our 16 gram gas cartridges.



Freeze limiter



Body



Phone ► +46 31-400 500
Fax ► +46 300-56 94 99

E-mail ► order@cryoiq.com
Website ► www.cryoiq.com



Gas cartridges with OxiBlock® technology

All our gas cartridges are equipped with a valve. The 16 and 25 gram gas cartridges fit CryoIQ® models DERM, DERM plus and PRO, and Cryoalfa® model SUPER and LUX. The gas cartridges have a built-in 5 µm filter and are filled with medical grade N₂O gas.

Order no.	Product description	Pack. size	Min order qty	Picture	Discount group
Gas cartridges					
CIQ-G-VN16	Gas cartridge 16 gram N ₂ O with valve for model CryoIQ® DERM, DERM plus, PRO and Cryoalfa® SUPER.	1	10 single package	G1	B
CIQ-G-VN25	Gas cartridge 25 gram N ₂ O with valve for model CryoIQ® DERM, DERM plus, PRO and Cryoalfa® SUPER.	1	10 single package	G1	B
CIQ-G-VN16x4	Gas cartridge 16 gram N ₂ O with valve for model CryoIQ® DERM, DERM plus, PRO and Cryoalfa® SUPER.	4	5 (4-pack)	G1	B
CIQ-G-VN25-LUX	New! Gas cartridge 25 gram N ₂ O with valve for model Cryoalfa® LUX.	1	10 single package		B



G1 / With valve
16/25 gram version



CIQ-G-VN16x4

CIQ-G-VN25



CIQ-G-VN16

OxiBlock® New!

The new OxiBlock® coating on the inside of the cartridge gives you a contamination-free environment for the gas. With this technique, we are the first brand that can offer you a cartridge that meets the new upcoming Medical Device Directive. What we have seen is that without our coating it is not possible to obtain a contamination free surface inside. This effects the durability and purity of the gas.

Warranty PLUS New!

With the meticulous control of all production steps of the gas cartridges, we are able to offer you, that if any faults in the devices from CryoIQ® or Cryoalfa® appears caused by use of CryoIQ® gas cartridges we will replace the unit with a new equal CryoIQ® device. Just remember it is important that you follow the cleaning instructions in our User Guide.

Filling New!

All of our cartridges are filled in-house. We have put a lot of efforts in developing a unique gas filling station that has the latest technique to optimize filling of cartridges. To establish a more precise and clean filling process, the filling station itself is equipped with nine separate 5 micron gas filters to ensure 100% purity. With our in-house filling process, we have 100% traceability of the gas used in our cartridges.

It was important for us to find an independent solution where we could fill and refill our cartridges. We have a close cooperation with Praxair who is our consultants and manufacturer of medical gas, but more importantly, we have the in-house knowledge. Today we have in our production two filling stations and are now looking for a fully automated filling line, with a capacity of more than 500 gas cartridges per day to meet future demands.

Medical gas New!

With our filling process, we can ensure that the medical gas we use is kept 100% clean during the entire filling process. We have seen not only once but several times that medical gas is filled in cartridges under conditions that does not fulfil medical standards, which renders the usage of medical gas pointless.

Cap New!

The fine and precise thread in the valve requires a precise cap to seal completely. Our new cap construction is in fine metal and it has its own sealing to ensure an even better shelf life for the cartridges. The previous solution was a combination of a plastic cap with a small plastic stopper that was placed inside the valve. In some cases we noticed that the cap did not seal properly. It also happened that the inner plastic stopper was forgotten to be removed when installing the device, which caused damages on the devices and a lot of frustration.


Shelf Life New!

Together with our gas manufacturer, Praxair, we have implemented a Just-In-Time system so that we can always offer you as our dealer the best shelf life for your gas cartridges. From the production day of the gas the durability is 3 years and with our system we never have gas older than 1 month in our production.

What we have seen before is that larger batches of production results in shorter shelf time for the gas cartridges. With our Just-In-Time system we have adjusted our production so that we never needs to build unnecessary stock. This system is also very cost effective which enables us to keep our prices on a low level.



If you need any specific picture or would like to have designed picture. Ask us and we will see how we can support you.



My notes

[illegible]




Phone ► +46 31-400 500
Fax ► +46 300-56 94 99

E-mail ► order@cryoiq.com
Website ► www.cryoiq.com



Marketing material

We offer following marketing material in smaller volume, for large quantities please contact us.
In our graphic studio we can produce custom designed material.



General



Aesthetic



Pediatric




Podiatric




Veterinary

Description: Leaflets for CryoIQ products
Format: Printed
Languages: English, German, French, USA
How to order: CIQ-F-MED-GEN-EN (DE/FR/US)
CIQ-F-MED-AES-EN (DE/FR/US)
CIQ-F-MED-PED-EN (DE/FR/US)
CIQ-F-MED-POD-EN (DE/FR/US)
CIQ-F-MED-VET-EN (DE/FR/US)
Option: Can be printed with Partner logo
Cost: Depending on volume

Description: Leaflets for CryoIQ products
Format: Digital - PDF in A4 size (4 pages)
How to order: Download from <http://cryoiq.com/downloads/>
Languages: English, German, French, Bulgarian, Russian, Spanish, Traditional Chines, more languages on the way.




Userguide medical



Userguide veterinary

Description: Userguides for CryoIQ medical products
Format: Printed
How to order: B-UG-MED-EN
Languages: Bulgarian, Danish, Dutch, English, Estonian, Farsi, Finnish, French, German, Greek, Hrvatska, Italian, Polish, Russian, Spanish, Swedish, Traditional Chines, more languages on the way.



Description: Rollup showing product portfolio
Format: RollUp 80 x 200 cm
How to order: CIQ-PRINT-ROLLPORT2016
Languages: Universal
Cost: €40 (approximately)



Description: DVD with pictures and other marketing material
Format: DVD
How to order: CIQ-DVD-2020
Languages: -
Cost: -